



# Changing Veterinary Inhalation Anaesthesia

**Optimum carbon dioxide absorption.**

**Irreversible colour change, when fully exhausted.** (Know when to change)

**Non-corrosive to skin.** (Safer to handle)

**Low resistance to gas flow.** (Isoflurane, Sevo)

**Low dust levels with consistent granule size.**

**Does not generate Carbon Monoxide CO under any clinical conditions.**

**Does not generate Compound A (CF<sub>2</sub>=(CF<sub>3</sub>)OCH<sub>2</sub>F) or any like compounds when used with Sevoflurane.**

**Does not generate Formic Acid when used with Sevoflurane.**



## ORDER INFORMATION:

AVP  
CDMV  
MVP  
VP  
WDDC

Amsorb Plus (Jug 4.5 kg)  
Amsorb Plus (Jug 4.5 kg)  
Amsorb Plus (Jug 4.5 kg)  
Amsorb Plus (Jug 4.5 kg)  
Amsorb Plus (Jug 4.5 kg)

Order# 9190102  
Order# 102180  
Order# 4020057  
Order# 1260018  
Order# 110533



**benson**  
medical industries inc.  
TORONTO • OTTAWA • CALGARY • VANCOUVER

151 esna park drive, unit 19  
markham, ontario L3R 3B1  
phone: (905) 475-0401  
fax: (905) 475-3656  
www.bensonmedical.ca

