





C3400 **MXR-Flowmeter
with Porter Scavenger Breathing Circuit**
3000 MXR-Flowmeter
4100 Three Litre Bag
5155-1 Porter Scavenger Breathing Circuit
• **C3400AV substitutes 5155-3AV Breathing Circuit**
See page 8 (with AVS)




C3420 **MXR Package C3400 plus**
2020 Solid Wall Mount
 -Telescoping
• **C3420AV substitutes 5155-3AV Breathing Circuit**
See page 8 (with AVS)

C3435-1 **MXR Package C3400 plus**
2035-1 Cabinet Mount
 -Slide Away
• **C3435-1AV substitutes 5155-3AV Breathing Circuit**
See page 8 (with AVS)




C3436 **MXR Package C3400 plus**
2036-3 Cabinet Mount
 -Slide
• **C3436AV substitutes 5155-3AV Breathing Circuit**
See page 8 (with AVS)

C3440 **MXR Package C3400 plus**
2040 Compact Mobile Stand
• **C3440AV substitutes 5155-3AV Breathing Circuit**
See page 8 (with AVS)



C3445-1 **MXR Package C3400 plus**
2045-1 Tall Mobile Stand
 for 2 + 2 "E" Tanks
• **3445-1AV substitutes 5155-3AV Breathing Circuit**
See page 8 (with AVS)

C3445E-1 **Same As Above plus**
1506-1 Positive Pressure/Demand Valve
• **C3445-1AV substitutes 5155-3AV Breathing Circuit**
See page 8 (with AVS)

 **NOTE:** When Scavenger is connected to a 6256-3 outlet with vacuum, order a 5602 found on page 5.
Alternative vacuum packages can be found on page 24. Necessary Gas Supply Tubing can be found on page 9.



CSM 2000 COMPACT
37 1/2" H x BASE 14" W x 16D



CSM 3000 PORTER CABINET MODULE
32"H x BASE 16" W x 18" D



CSM 4000 TOWER
51" H x BASE 21" W x 18"D

Our latest flush mount nitrous oxide sedation mixer, incorporating all of the safety and design features engineered into our popular PORTER MXR Vertical Flowmeter is now available in mobile cabinets aesthetically designed to meet the needs of the modern dental operator.

- CSM 2000

Porter Compact Portable
Accommodates 2 "E" oxygen cylinders and 1 "E" nitrous oxide cylinder.
Complete less demand resuscitator
- CSM 3000

Porter Cabinet Module
Nitrous oxide sedation system, accommodates 2 "E" oxygen cylinders and 1 "E" nitrous oxide cylinder. Complete, less demand resuscitator
Dental cabinet not included
- CSM 4000

Porter Tower Portable (offers pipeline economy)
Accommodates 1 "M" oxygen cylinder and 1 "E" nitrous oxide cylinder.
Complete less demand resuscitator
Includes:
Porter MXR-1 Cabinet Mounted Mixer
Porter Autoclavable Gas Scavenging Rubber Goods with 5054A adult hood assembly
Gas Scavenging Interface with adaptor for H.V.C.
Two-Stage Regulators for O₂ and N₂O
Deluxe Cabinet, arborite finish with carpet casters
- 1506-1

CPR/Demand Valve Resuscitator with 1575 Adult Mask
6' hose with Porter A-1675 quick connect
- 1575-1

Universal Size Full Face, Clear Mask
- 5054A

Porter Adult Hood Assembly with 3 inserts
- 5454B

Pedo Hood Assembly with 3 inserts

Economy of use and purchase savings can be achieved by selecting a portable of choice as an option to an expensive medical gas pipeline system.