

## CPR/Demand Valve Resuscitator

The PORTER/O-TWO Demand Valve Resuscitator provides a fast, simple and effective means to ventilate non-breathing patients and provides 100% oxygen "on demand".

Designed to meet the varying requirements of the Emergency Medical Services/Rescue Personnel, it is the only Demand Valve Resuscitator that comes standard with an Audible Pressure Relief System and sensitive Manual Ventilation Button, conveniently located on the front of the resuscitator. The Manual Ventilation Button helps to ensure a patient airway by enabling both hands to remain on the mask to facilitate a superior mask to face seal (without having to deal with cumbersome trigger levers or inconveniently located top buttons).

The Manual Ventilation Flow rate is fixed at 40 lpm in line with the 1992 American Heart Association "Recommendations and Guidelines for Resuscitation."

During spontaneous respirations, the Demand Valve Resuscitator is designed to provide 100% oxygen to a breathing patient with minimal respiratory effort. The oxygen demand flow is controlled by the patient's inspiratory effort, with flow rates up to 160 lpm, to accommodate their own respiratory distress situation.



- |        |  |
|--------|--|
| 1506-1 | <b>CPR/Demand Valve</b> , only<br>with 1575 adult mask, 6 ft. hose,<br>and PORTER A-1675 quick connect   |
| 1536-1 | <b>CPR/Demand Valve</b> , mobile cart<br>with 1575 adult mask, 10 ft. hose, regulator and yoke<br>(excludes cylinder) and PORTER A-1675 quick connect                                |
| 1540-1 | <b>CPR/Demand Valve</b><br>in Softpack Clamshell Carrying Case, including<br>flow control regulator and wrench ( <i>excludes cylinder</i> )  |
| 1542   | <b>Oxygen Therapy Inhalator Kit</b><br>in Softpack Clamshell Carrying Case with adjustable<br>1-25 lpm regulator, disposable oxygen mask with<br>wrench ( <i>excludes cylinder</i> ) |
| 1544   | <b>Bag-Valve-Mask Resuscitator</b><br>Silicone, complete with bag, valve, universal mask,<br>4 airways, reservoir system and carrying case   |
| 1575-1 | <b>Universal Size Full Face, Clear Mask</b>  |
| A-1675 | <b>Porter O<sub>2</sub> diss/Quick Connect</b>   |



1506-1



1544



Other options for oxygen administration equipment are available.  
Please consult your Benson Medical Representative to discuss.