



8" Automatic Steam Sterilizer Specifications:

Chamber size:	8 "
Overall Dimensions:	18.5" W x 20" D x 14.5" H
Electrical:	Supply: 117V +/- 10%, 60 Hz
	Loading: 117V, 1400w
	Fuses: 117V, 15A, 800mA
Weight:	Net: 64 lbs.
	Shipping: 70 lbs.
SES-2000E	Automatic Steam Sterilizer System with 3 trays, 8" chamber
SES-1300	Tray Package includes three (3) small trays and one (1) lifter
SES-1200	Pouch Rack for 8" chamber

Note: Porter Instrument reserves the right to make changes and improvements to this equipment without notification.



12" Automatic Steam Sterilizer Specifications:

Chamber size:	12"
Overall Dimensions:	23.6" W x 24.5" D x 20.6" H
Electrical:	Supply: 240V +/- 5%, 60 Hz or 208V +/-5%
	Loading: 240V, 3Kw
	Fuses: 250V, 0.2A
Weight:	Net: 153 lbs.
	Shipping: 173 lbs.
Optional Printer:	Requires thermal type printing paper, 58mm wide
SES-2012E	Automatic Steam Sterilizer System with 5 trays, 12" chamber
SES-2012EP	Automatic Steam Sterilizer System with Optional Printer, 5 trays, 12" chamber
SES-1212	Pouch Rack for 12" chamber

Note: Porter Instrument reserves the right to make changes and improvements to this equipment without notification.