



**NOTE: Modular Cable should NOT be cut or spliced.  
The voltage drop is too great and it will affect the system operation.**

STOCK NO.	DESCRIPTION
A. A-2554-000	<b>O2 Gauge, 0-150 psi, 2" dial</b>
B. A-2554-001	<b>N2O Gauge, 0-150 psi, 2" dial</b>
C. B-1322	<b>O2 diss Check Valve Assembly</b>
D. B-1323	<b>N2O diss Check Valve Assembly</b>
E. B-1324	<b>Pressure Relief Valve Assembly</b>
F. C-394-000	<b>Pressure Switch Assembly</b>
G. A-2144	<b>diss/Copper Pigtails (straight)</b> <b>Oxygen - 1/2" O.D. copper x diss O<sub>2</sub> female</b> <b>Nitrous Oxide - 3/8" O.D. copper x diss N<sub>2</sub>O female</b>
H. A-2041	<b>Sentinel Interconnect Board Assembly</b>
I. D-160-PL	<b>Molded Plastic 3222MX Cover</b>
1. 5775	<b>75' Sentinel Modular Connector Cable (standard)</b>
5700-150	<b>150' Sentinel Modular Connector Cable (special order)</b>
5700-200	<b>200' Sentinel Modular Connector Cable (special order)</b>
2. 5700	<b>50' Sentinel Modular Connector Cable</b>
3. 5700-1T	<b>100' Teflon Sentinel Modular Connector Cable</b>
4. KC13161-000	<b>Power Supply for Wall Sentinel</b>
6. C-637-000	<b>Power Supply for Desk Sentinel</b>
7. A-2083	<b>Modular Cable Connector Box</b>
PA-636	<b>Sentinel Light Bulb-" II Glass"</b>
PA-571	<b>Sentinel Light Bulb-" Id Bayonet Style"</b>