



Items Below in GRAY Boxes are STEAM AUTOCLAVEABLE!

STOCK NO.	DESCRIPTION
1. 5054A	Porter Adult Hood-Complete with 3 Liners
1a. 5054-1	Adult Inner Hood Liners-Package of 3
1b. 5054C	Porter Adult Hood-with 1 Liner
2. 5054B	Porter Pedo Hood-Complete with 3 Liners
2a. 5054-2	Pedo Inner Hood Liners-Package of 3
2b. 5054D	Porter Pedo Hood-with 1 Liner
3. 5056	Coaxial Inhaler Tubing
4. 5057	Tension Slide Adjuster
5. 5058	Fresh Gas "Y" Connector
6. 5059	Vacuum Hose (8 feet)
PA-679-000	Vacuum Hose by the foot
7. 5060A	Fresh Gas Corrugated Tubing (3 feet)
5060B	Fresh Gas Corrugated Tubing (6 feet)
8. 1570-22	22mm Right Angle Adapter
1570-24	24mm Right Angle Adapter
9. 5061	Mask to Tubing Plastic Connectors (set of 2)
10. 5062	Vacuum Line "Y" Connector
11. 5501-RK	In-Line Vacuum Kit (can be utilized as pictured above)
12. 4100	3 Litre Bag
4100-3LF	3 Litre Bag Latex Free
4100-2LF	2 Litre Bag Latex Free
13. 5063	"T" Adapter for In-Line Vacuum Block
14. 5064	"Straight" Adapter for In-Line Vacuum Block
15. 5065	Vacuum Tube Holder
16. 5066	Porter Liner Valve Only
17. 5054K	Porter Liner Connector Kit (pkg of 6)-for liners WITH removable connectors ONLY (old style)
5602	Vacuum Quick Connect Disposable Liners found on page 8.