






AVS-5000	<b>Automatic Vacuum Switch System</b>	
5501	<b>Wall or Cabinet Mount</b> for Hook-Up to HVE Line includes:	
5450A	Vacuum block w/flow tube, wall, cabinet post or mobile stand mount	
5320	Wall or cabinet bracket	
A-1911	10 foot 3/8" Vacuum tubing w/fittings	
A-1907	HVE adapter for 3/8" tubing	
5651	Right angle adapter	
5521	<b>Outlet Station Mount</b> includes:	
5450C	Vacuum block w/flow tube, outlet station mount	
5651	Right angle adapter	
5450A	<b>Vacuum Block Assembly</b> wall or post mount	
5450A-RK	<b>Vacuum Block Retrofit Kit</b> converts 5450A to 5450C	
5450C	<b>Vacuum Block Assembly</b> outlet station mount	
5450B-RK	<b>Vacuum Block Retrofit Kit</b> converts 5450B to 5450C	
A-2095-BR	<b>Flow Nozzle</b> (Brass Filter/Restrictor)	
A-1902	<b>Acrylic Flow Tube Assy</b> for Vacuum Block	
A-1907	<b>HVE Adapter Assy.</b> for 3/8" polyflow tubing includes: Polyflow Ferrule and Nut	
A-1911	<b>10 foot 3/8" Vacuum Tubing</b> with Fittings	
5651	<b>Right Angle Adapter</b>	
PA-296-003	<b>Viton Sealing Washer for 5651</b>	
5320	<b>Wall or Cabinet Bracket</b> Holder for 5450A	

## Oxygen Resuscitation

A-1675	<b>Porter O<sub>2</sub> diss/Quick Connect</b>	
1560	<b>Regulator and Yoke</b> w/diss Male Connection	
1566	<b>Regulator and Yoke</b> w/Porter Imperial Connection	
A-1649-000	<b>"T" Handle for 1560 and 1566</b>	