

**2450WF MXR-1 Flushmount Flowmeter**  
(Bottom Connected)



**2455WF MXR-1 Flushmount Flowmeter**  
(Back Connected)

**2465WF MXR-1 Flushmount Flowmeter**  
(Universal connection)(back)  
One model with connections that are easily switched in the field from back to bottom

**2470WF same as above (bottom)**

- Packages above include flowmeter with bag tee, Porter scavenger breathing circuit, vacuum switch and complete flexible installation kit.
- Flexible installation is the more prevalent type. Hard copper installation and bezel mounting bracket are required when flowmeter is mounted in a wall. When hard copper installation kit is required, please specify by ordering 2450WH.

## Flexible Installation



- All components are DISS (Diameter Index Safety System) and interchangeable to accommodate the maximum number of installation requirements.