



8065* **Hose, O₂/N₂O, diss/Porter, 5 ft.**




8055* **Hose, O₂, diss/Porter, 5 ft.**



8005* **Hose, O₂, diss/diss, 5 ft.**




8505* **Hose, N₂O, diss/diss, 5 ft.**



8015* **Hose, O₂, diss/NCG (Chemetron), 5 ft.**

8515* **Hose, N₂O, diss/NCG (Chemetron), 5 ft.**



8025* **Hose, O₂, diss/Puritan, 5 ft.**


8525* **Hose, N₂O, diss/Puritan, 5 ft.**

for use with McKesson & Fraser Harlake
outlet stations




8035* **Hose, O₂, diss/Schrader, 5 ft.**

8535* **Hose, N₂O, diss/Schrader, 5 ft.**



8045* **Hose, O₂, diss/Oxequip, 5 ft.**

8545* **Hose, N₂O, diss/Oxequip, 5 ft.**



8075* **Hose, O₂, diss/Ohio Medical/ Matrx, 5 ft.**

8575* **Hose, N₂O, diss/Ohio Medical/ Matrx, 5 ft.**

* ALL STOCK NUMBERS LISTED ABOVE ARE FOR A 5 FOOT
HOSE ASSEMBLY. Last digit indicates length in feet.
Example: 8053 - Hose, O₂, diss/Porter - 3 feet long.

**Note: Hoses longer than 5 feet are custom made.
ALL SALES ARE FINAL.**