



BLACK GUEDEL AIRWAYS



20742

20736

20700



FLEXCONNECTOR

Guedel Airway

The Guedel airway provides a large tubular lumen for increased air exchange during insertion of suction apparatus. These single use items are available in black rubber, clear polyvinyl and semi clear polyethylene plastic.

SIZE	CLEAR*	BLACK RUBBER
000, 30mm	11082	11066
00, 40mm	11083	11067
0, 50mm	11084	11069
1, 60mm	11085	11071
2, 70mm	11086	11073
3, 80mm	11087	11075
4, 90mm	11088	11077
5, 100mm	11089	11078
6, 110mm	11090	11080

SIZE	SEMI CLEAR*	COLOUR CODED*
Infant, 50mm	11700	11720 BLUE
Small Child, 60mm	11702	11722 BLACK
Child, 70mm	11704	11724 WHITE
Small Adult, 80mm	11706	11726 GREEN
Medium Adult, 90mm	11708	11728 YELLOW
Large Adult, 100mm	11710	11730 RED

10 per package *Latex-free

Flex Connector Components

MALE ADAPTOR 15mm male	DELRIN	20734
	STAINLESS STEEL	20736
FEMALE ADAPTOR 15mm female	DELRIN	20740
	STAINLESS STEEL	20742

Semi-Clear Lightweight

Flexible 9/16" smooth I.D. tubing with 15mm male and female ends.

6" overall length	20700
9" overall length	20702
12" overall length	20704

King Flexconnectors **DISPOSABLE**

Accordian-style tubing provides flexibility in the placement of the anaesthesia circuit during difficult procedures such as facial or head surgery.

- Contracts or expands to suit requirements
- Clear corrugated tubing improves visualization
- Individually wrapped in case of 50.
- 15mm Female x 15mm Male ends.

4" to 6 1/2" length	101035
6 1/2" to 14" length	1010412