




10043



20804



11673



10010



ANTI-STATIC TUBING



10131



10141



10135



KAB-001

## THREE PORT DELRIN ADAPTOR ..... 10043

Has 19mm female port and two 19mm male ports.

## SMALL CYLINDER WRENCH ..... 20804

## BODAK SEAL (Not Shown) ..... 10020

## TEST LUNG ..... 10010

Durable 500 cc test lung with 22mm female bushing.

## TUBE TREE ..... 11673

A versatile breathing tube support which accepts adult and pediatric corrugated breathing tubes. It may be placed either above or below the table and supports tubes whether patient is in the supine, prone or lateral position. Prevents drag and pull on mask or E/T tube and helps prevent accidental extubation, advancement of E/T tube to an endobronchial position, or kinking of the E/T tube at the connector end.

## ANTI-STATIC TUBING

Manufactured from high quality black latex rubber. Anti-static with yellow stripe.

## TUBING 4.75 x 1.6mm ..... 10131

(I.D. 3/16" x Hose Thickness 1/16") Roll - 20m

## TUBING 4.75 x 3.2mm (Not shown) ..... 10133

(I.D. 3/16" x Hose Thickness 1/8") Roll - 20m

## TUBING 6.35 x 1.6mm ..... 10135

(I.D. 1/4" x Hose Thickness 1/16") Roll - 20m

## TUBING 6.35 x 3.2mm (Not shown) ..... 10139

(I.D. 1/4" x Hose Thickness 1/8") Roll - 20m

## TUBING 5.5 x 1.6mm ..... 10141

(I.D. 7/32" x Hose Thickness 1/16") Roll - 20m

## TUBING 5.5 x 3.2mm (Not shown) ..... 10145

(I.D. 7/32" x Hose Thickness 1/8") Roll - 20m

## King Systems KAB Disposable CO<sub>2</sub> Absorber

- Ideal for potentially infectious procedures
- Protects both patient and equipment from cross contamination
- Totally self-contained. Includes clear housing for visualization of carbon dioxide absorbent, APL valve, unidirectional valves and bag mount
- Easily mounted on all models of anaesthesia machines
- Latex-free materials

## KING SYSTEMS KAB CO<sub>2</sub> ABSORBER ..... KAB-001