




VS8QSENEZ



VAPORIZER POLE MOUNT



12020


Penlon Sigma Delta Inhalation Vaporizer

The Penlon Model Sigma Delta Inhalation Vaporizer delivers superb performance under varying flowrates and temperature.

KEY FEATURES

- Ten Year Service Free performance (Halothane @ five year recommended service intervals)
- Available mounting - Drager, Narkomed, Cagemount and Selectatec®
- Accurate delivery of vaporizer concentrations, particularly at low flows
- Approximate weight - 5 kg
- Holds approximately 250 ml of anaesthetic agent

**Variation of output with flow rate (Temperature: 20°C)
Isoflurane**



MOUNTING ARRANGEMENT

	CAGEMOUNT 23MM TAPER	SELECTATEC	DRAGER JULIAN PLUG-IN	NARKOMED
Sevoflurane 0 to 8% vol	VS8QCENEZ	VS8QSENEZ	VS8QDENEZ	VS8QNEZ
Isoflurane 0 to 5% vol	VI5KCENEZ	VI5KSENEZ	VI5KDENEZ	VI5KNEZ
Halothane 0 to 5%vol	VH5KCENEZ	VH5KSENEZ	VH5KDENEZ	VH5KNEZ
Enflurane 0 to 7% vol	VE7KCENEZ	VE7KSENEZ	VE7KDENEZ	VE7KNEZ

Note: All of the above vaporizers are equipped with either Key Fill adaptors or Quik Fil® adapter (Sevoflurane only)

VAPORIZER POLE MOUNT BMVAPMT

VAPOFIL™ Keyed Anaesthetic Agent Adaptor

Available for use with: Isoflurane • Enflurane • Halothane

Vapofil™ Keyed Filler - Isoflurane (Purple).....	12016
Vapofil™ Keyed Filler - Enflurane (Orange)	12018
Vapofil™ Keyed Filler - Halothane (Red)	12020