

Stethoscopes and Accessories



Stethoscopes

BINAURAL STETHOSCOPE 11120

Anaesthesia length-tubing from "Y" to female Luer connector is 36" long to allow for movement around patient.

BINAURAL STETHOSCOPE 11122

Same as #11120 but with Luer-Lock fitting

Chestpieces

WENGER

Heavyweight chest piece of chrome plated brass. May be used on some areas of the body without being taped. Includes male Luer fitting for attachment.

Adult **11060**

Child **11061**

MARTIN 11059

Lightweight Martin type accommodates standard plastic I.V. tubing. The angle of the inlet keeps the chest piece flat.

INFANT SIZE 11057

Lightweight, with male Luer fitting.

PREMATURE SIZE 11058

Lightweight, with male Luer fitting.

Accessories

STANDARD BULB AND VALVE 22208

Inflation bulbs for use with blood pressure cuffs and pocket aneroids. Includes moulded latex bulb and metal valve.

BLOOD PRESSURE INFLATION VALVE ONLY 22210

BLOOD PRESSURE INFLATION BULB ONLY 22212

CONNECTOR SET 22214

Male and female Luer.

CONNECTOR SET 22216

Male and female Luer lock.

BLOOD PRESSURE TUBING (Not Shown) 22218

Non-conductive 3/16" I.D. x 3/32" Wall.

ONE HAND ANEROID MANOMETER 22209

Includes manometer, valve and bulb

ANEROID BLOOD PRESSURE GAUGE ONLY 22206

0 to 300mmHg Gauge