

I.V. Pole Stands

KEY FEATURES:

- Patented GFF™ Non-Stain Base resists most chemicals found in hospitals.
- Squeeze Release Latch with backup safety lock.
- Duo-Casters that are easy to clean.
- Four-layer nickel and chrome plating on all steel

SPECIFICATIONS:

- Base Width: 24"
- Maximum Height: 100"
- Minimum Height: 52.5"
- Caster Diameter: 2"
- Weight: 14 lbs.
- Pole Diameter: 1"
- Five leg stand

I.V. Chrome Stat-Stand™ with 2-Prong Hook **1050-2**

I.V. Chrome Stat-Stand™ with 4 Prong Hook **1050-4**

Combination Rolling Stand with Tray

- 10 1/2" x 12" steel tray
- Two sets of strap adjustment slots
- Vel-cro® Hold-down strap, 1" diameter, 4 feet long
- Cord Wrap Bracket
- Guidance handle
- Five leg stand
- 39" overall height
- Pole diameter 1 3/4"

Rolling Stand with tray **1064**

ACCESSORIES:

Basket with Pole Mount (12" long x 6" wide x 6" deep) **1034**

Basket with Pole Mount (12" long x 12" wide x 8" deep) **1035**

Guidance Handle (Fits 3/4" to 1 1/4" diameter poles) **1048**

Dual Pump Holder (Uprights are 1" dia. X 12" high) **1049**

