



3601

Sechrist Air Oxygen Mixer

The Sechrist Model 3500HL High/Low Flow Air Oxygen Mixer is a proven precision proportioning device designed to accurately mix and deliver oxygen and medical air in selected FIO2 concentrations from .21 to 1.0.

KEY FEATURES:

- Years of proven performance
- Balance regulator design allows accurate delivery of gases even if supply gas pressure fluctuates.
- Includes a 0.1 micron air filter as a standard accessory on all Sechrist Mixers
- Simple to use and maintain
- Accurate to 3% of scale
- Bleed flow of 8 lpm
- Dual Oxygen DISS outlets
- Available in both wall and pole mount configurations
- Bypass alarm system should either supply gas fail

Sechrist Air Oxygen Mixer w/ pole mount	3600
Sechrist Air Oxygen Mixer with rear wall mount.....	3601



15408-03

Timeter® Soft-Touch Oxygen Flowmeter

- Durable, Dual tapered metering tube for outstanding accuracy
- Five-year warranty
- Delivers oxygen flows from 0 to 15 lpm, and flush flows of over 70 lpm
- Backpressure compensated

Timeter O2 Flowmeter	15408-03
Includes DISS Female Inlet Adaptor with Hand Tight Nut	

Oxygen Regulators

- Chrome-plated forged brass construction
- Safety relief valve
- Preset 50 psi outlet pressure
- Cylinder contents gauge



50600



50601

Oxygen Regulators with CGA540 Oxygen Nut and Nipple.....	50600
Oxygen Regulators with CGA870 Yoke Assembly	50601

O-Two Systems Disposable Easy Grip® Resuscitator

- Rapid bag recoil feature facilitates high ventilation rates
- Standard 15/22mm swivel patient connector.
- Standard pressure limiting device on child and infant models
- Available in cases of 12

EASY GRIP® DISPOSABLE RESUSCITATOR

Adult, w/ mask, reservoir system and O2 tubing	01BM3500-Cs
Child, w/ mask, reservoir system and O2 tubing	01BM3510-Cs
Infant, w/ mask, reservoir system and O2 tubing	01BM3520-Cs

NEW!

O-Two Systems Resuscitation Smart Bag®

The SMART Bag® (Synchronized Manual Actuation Response Technology) has a unique actuating mechanism inside the neck bushing of the bag that limits the excessive flow of gas from the bag into the patient's airway, significantly reducing the risk of barotrauma, gastric insufflation and subsequent complications.

KEY FEATURES:

- Visual indicator reminds the operator to reduce force being applied to the bag
- Compensates for patients with poor compliance by allowing higher delivered pressures
- All bag/valve/masks are latex-free
- Child Bag/Valve includes pressure safety relief valve
- Minimizes risk of gastric insufflation and barotraumas
- Attachment port for supplemental oxygen

Reusable Silicone SMART BAG®, Adult with Facemask and Reservoir System	01BM3100
Reusable Silicone SMART BAG®, Child with Facemask and Reservoir System.....	01BM3110
Disposable SMART BAG®, Adult w/ Facemask & Reservoir System (6/CASE).....	01BM3200-Cs
Disposable SMART BAG®, Child w/ Facemask & Reservoir System (6/CASE).....	01BM3210-Cs

Mobile Medical Gas Cylinder Carts

Fits all "D" or "E" size medical gas cylinders.

Mobile Cylinder Cart, holds 12 cylinders (FOUR WHEELS)	MC-12
Mobile Cylinder Cart, holds 6 cylinders (FOUR WHEELS)	MC-6
Mobile Cylinder Cart, holds 2 cylinders (TWO WHEELS)	MC-02
Mobile Cylinder Cart, holds 1 cylinder (TWO WHEELS)	MC-01



EASY GRIP



SMART BAG



MC-02



MC-01