


## Feline Orogastric Tube



### DESCRIPTION

- A closed-end single lumen silicone catheter with distal sideholes used for the administration of enteral feeding solutions.
- Proximal fitting accepts a tapered syringe adaptor.
- Recommended for VETERINARY USE ONLY.


### ORDER No.

### Remarks

V-FOT-18-43

Silicone Orogastric Tube, 18 French 43cm length. Autoclavable.

## Feline Oesphagostomy Tubes



### DESCRIPTION

- A closed-end single lumen silicone catheter with distal sideholes placed surgically for the administration of enteral feeding solutions.
- Removable proximal fitting accepts standard Luer lock syringes.
- Recommended for VETERINARY USE ONLY.

### ORDER No.

### Remarks

V-FET-14-23

A 14 French silicone oesphagostomy tube, 23cm length with removable proximal fitting. Catheter only is autoclavable.

V-FET-19-38

A 19 French silicone oesphagostomy tube, 38cm length with removable proximal fitting. Catheter only is autoclavable.

V-PFLLA-VTA

Female Luer lock tapered adaptor only.

V-PFLLP-100

Luer lock plug only.

## Feline Mouth Gag

### ORDER No.

### Remarks

V-FMG-100

A reusable mouth gag designed to assist in the introduction of orogastric tubes, up to 18 French outside diameter.



### REFERENCE

Dr Gary Norsworthy, DVM, DABVP, San Antonio, Texas, USA.