


# Drainage

## Chest Drainage Tubes – Silicone



### DESCRIPTION


- Open-end single lumen silicone catheter with large multiple sideholes over the distal portion of the catheter.
- The catheter can be used for the drainage of air or fluid from body cavities including the pleural space of the thoracic cage.
- Markings at 2.0 and 4.0cm from the proximal sidehole assist in correct catheter placement.
- Intended for VETERINARY USE ONLY.



ORDER No.	Remarks
V-PCT-12	12 French straight drainage tube with 2.2mm I.D. and 4.0mm O.D. Sideholes over distal 5.0cm of catheter.
V-PCT-14	14 French straight drainage tube with 2.5mm I.D. and 4.7mm O.D. Sideholes over distal 7.0cm of catheter.
V-PCT-20	20 French straight drainage tube with 4.0mm I.D. and 6.6mm O.D. Sideholes over distal 10cm of catheter.
V-PCT-30	30 French straight drainage tube with 7.0mm I.D. and 10.0mm O.D. Sideholes over distal 12cm of catheter.

**Note:** V-PFLLA-VTA (page 46) is recommended for use with V-PCT-12, V-PCT-14 and V-PCT-20  
 V-PFLLA-VTA-L (page 46) is recommended for use with V-PCT-30

## Metal Stiffening Stylets



Stainless Steel Stiffening Stylet



ORDER No.	Remarks
V-TCS-52-40	Metal stiffening stylet, 1.3mm O.D., 40cm long. Recommended for 12 French drainage tubes.
V-TCS-90-40	Metal stiffening stylet, 2.3mm O.D., 40cm long. Recommended for 14 French drainage tubes.
V-TCS-119-40	Metal stiffening stylet, 3.2mm O.D., 40cm long. Recommended for 20 French drainage tubes.
V-TCS-187-45	Metal stiffening stylet, 4.7mm O.D., 45cm long. Recommended for 30 French drainage tubes.