



Pneumothorax Sets



Catheter Support Fixation Disc (supplied separately)



Catheter and Needle

DESCRIPTION

- Used to treat a spontaneous or tension pneumothorax.
- The Luer lock fitting on the catheter can be connected to suction apparatus or a Heimlich valve.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

Canine

ORDER No.

Remarks

V-TPT-1020U

10 French, 40cm long polyurethane catheter with 20cm long introduction needle stylet.

V-FD-100

Catheter support low profile fixation disc.

Feline

ORDER No.

Remarks


V-TPT-815U

8.3 French, 35cm long polyurethane catheter with 15cm long introduction needle stylet.

V-FD-100

Catheter support low profile fixation disc.

Accessory Items



Heimlich Chest Drainage Valve



Drainage Connecting Tube

ORDER No.

Remarks

V-CASP-A-FORD

Heimlich one-way valve and connecting tube.

V-CTU14.0-30

Connecting tube only.