



DESCRIPTION

- A 16 French chest drain suitable for use in medium to large dogs.
- The chest tube is introduced with the supplied trocar using a direct puncture technique.
- Multiple large sideports facilitate drainage of pus, haem or chyle.
- A radiopaque stripe allows tube position to be checked radiographically.
- The Luer lock fitting on the catheter can be connected to suction apparatus.
- Intended for VETERINARY USE ONLY.

ORDER No.

Remarks

V-TPTL-1630V

16 French, 30cm long catheter with 12cm long introduction needle stylet.

What Chest Drain Do I Need?

PNEUMOTHORAX SET

This catheter has been designed to treat spontaneous or tension pneumothorax in dogs or cats. It has multiple small sideholes suitable for drainage of air and fluids of low viscosity. It has been used to drain selected hydrothorax, haemothorax and chylothorax cases.

The pneumothorax catheter is introduced by direct puncture into the thoracic cavity, the catheter is advanced off the stylet into the desired position. This allows for rapid and accurate placement.

CHEST DRAINAGE TUBES — SILICONE

This catheter comes in a variety of sizes suitable for dogs and cats. It is introduced by blunt dissection and can be stiffened with a stylet to facilitate placement if desired. It is open-ended, has large sideholes and is suitable for drainage of viscous fluids such as pus. It is ideal for pyothorax in cats. Because it is silicone, it is reautoclavable and reusable. The catheter comes without a fitting. V-PFLLA-VTA multi-purpose tubing adaptor is suitable for V-PCT-12, V-PCT-14, V-PCT-20 tubes and should be ordered separately.

CHEST DRAINAGE TUBES WITH TROCAR

This chest drain is suitable for use in dogs and has large sideholes capable of draining viscous fluids such as pus. It is introduced by direct puncture into the thoracic cavity. Once the distal end is within the thoracic cavity, the catheter is advanced off the stylet into the desired position. This allows rapid accurate placement.

This chest drain provides an alternative to the silicone chest drains V-PCT-14 and V-PCT-20 for those veterinarians who prefer a direct puncture to a blunt dissection introduction technique. It is intended for one-time use only.