


"Grenade" Continuous Suction Bulb



100cc Bulb

150cc Bulb

400cc Bulb

DESCRIPTION

- Allows short-term suction to be applied to body cavities enabling the drainage of either air or fluid.
- Bulb is ideal for use in small animals.
- One-way valve prevents the return of air or fluid back into the drained cavity.
- Clear suction bulb allows easy visualisation of collected fluids.
- Graduations facilitate greater accuracy in the estimation of fluid loss.
- Luer lock allows the device to be securely attached to drainage catheters with proximal Luer lock fittings.
- Side tab allows the device to be anchored to the body by either bandage or tape.
- Supplied sterile in peel-open packaging.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.



ORDER NO

REMARKS

V-GSB-100	100cc suction bulb with one-way valve and Luer lock fitting.
V-GSB-150	150cc suction bulb with one-way valve and Luer lock fitting.
V-GSB-400	400cc suction bulb with one-way valve and Luer lock fitting.