



"Flexi-drain" Multi-Purpose Wound Drain – Silicone




DESCRIPTION

- A silicone tubular drain that may be trimmed or used as one piece for surgical wound drainage. Drain can be split along its length or trimmed to suit individual requirements.
- Catheter material is radiopaque.
- Intended for VETERINARY USE ONLY.

ORDER No.	Remarks
V-PWD-12-40	40cm long drain consisting of twelve (12) catheter lumens joined together, each of 3.0mm diameter.



Silicone Penrose Drain



DESCRIPTION

- A silicone penrose drain used to promote surgical wound drainage.
- Catheter is radiopaque to facilitate placement verification.
- Silicone softness may minimise wound irritation and trauma. It is suitable for use in latex-sensitive patients.
- Non-reactive material which can be autoclaved.
- Intended for VETERINARY USE ONLY.

ORDER No.	Remarks
V-SPD-003000	18mm x 6.0mm diameter, 30.5cm long drain.
V-SPD-002000	9mm x 6.0mm diameter, 20.3cm long drain.