



## DESCRIPTION

- A closed endhole silicone catheter with distal balloon cuff used for intermittent drainage of urine from the bladder or contrast studies of urogenital tract.
- The balloon aids in retaining the catheter inside the bladder.
- Wire guide stylets are available to assist in catheter introduction.
- Intended for VETERINARY USE ONLY.



**Please note:** Only sterile water to be used to inflate balloon

ORDER No.	Remarks
V-PFC5-30	5.0 French Foley, 30cm long (12") with 0.9mm I.D. and 1.7mm O.D.
V-PFC5-55	5.0 French Foley, 55cm long (22") with 0.9mm I.D. and 1.7mm O.D.
V-PFC6-30	6.0 French Foley, 30cm long (12") with 1.1mm I.D. and 2.0mm O.D.
V-PFC6-55	6.0 French Foley, 55cm long (22") with 1.1mm I.D. and 2.0mm O.D.
V-PFC8-30	8.0 French Foley, 30cm long (12") with 1.3mm I.D. and 2.7mm O.D.
V-PFC8-55	8.0 French Foley, 55cm long (22") with 1.3mm I.D. and 2.7mm O.D.
V-PFC10-30	10 French Foley, 30cm long (12") with 2.0mm I.D. and 3.3mm O.D.
V-PFC10-55	10 French Foley, 55cm long (22") with 2.0mm I.D. and 3.3mm O.D.

## WIRE GUIDE INTRODUCTION STYLETS

V-SF-25-60	Wire guide stylet for 5.0 and 6.0 French, 55cm (22") Foley catheter.
V-SF-38-60	Wire guide stylet for 8.0 French, 55cm (22") Foley catheter.
V-SF-45-60	Wire guide stylet for 10 French, 55cm (22") Foley catheter.