


Central Venous Catheter Set with Seldinger Introduction Technique



DESCRIPTION

- Used for venous pressure monitoring, blood sampling and administration of drugs and fluids.
- All catheters introduced using the Seldinger introduction technique. See previous page for instructions.
- All set components available upon request. Other sizes and lengths available upon request.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.



Canine

ORDER No.

Remarks

V-PUM-501J-V15

Set consists of 5 French polyurethane catheter, 15cm long, 6 French vessel dilator, 50cm "J" wire guide and 18 gauge, 7.0cm long needle.

Feline

ORDER No.

Remarks

V-PUM-401J-V10

Set consists of 4 French polyurethane catheter, 10cm long, 5 French vessel dilator, 50cm "J" wire guide and 20 gauge, 4.0cm long needle.

V-PUM-301J-V8

Set consists of 4 French polyurethane catheter, 8cm long, 30cm "J" wire guide and 21 gauge, 4.0cm long needle.