




Fluid Therapy

Endo-Sof™ Subcutaneous Catheter Set



DESCRIPTION

- An Endo-Sof™ catheter with multiple fenestrations used for the administration of subcutaneous fluids.
- Catheter introducer allows easy placement.
- The short stylet enhances the anchorage of the catheter as well as occluding the female Luer lock.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.



ORDER No.	Remarks
V-FSQ-7-12	7 French, 12cm long fenestrated subcutaneous catheter set.
V-FSQ-7-18	7 French, 18cm long fenestrated subcutaneous catheter set.
V-FSQ-7-25	7 French, 25cm long fenestrated subcutaneous catheter set.
V-FSQ-40-4-SW	Short stylet, ordered separately (shown above).