


Disposable Intraosseous Infusion Needles



DESCRIPTION

- Used as an alternative to intravenous catheter access during emergencies.
- Positioning of the needle within the tibial, femoral or ulnar medullary cavity facilitates infusion of resuscitative drugs and fluid.
- Supplied in peel-open package.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

ORDER No.	Remarks
V-DIN-14-3.0-POWCH	14 gauge needle, 3.0cm length.
V-DIN-16-3.0-POWCH	16 gauge needle, 3.0cm length.
V-DIN-18-3.0-POWCH	18 gauge needle, 3.0cm length.
V-DIN-20-3.0-POWCH	20 gauge needle, 3.0cm length.