


Biopsy Devices

Vet-Core Biopsy Needles



DESCRIPTION

- Used for soft tissue core biopsy.
- Light-weight and compact, Vet-Core Biopsy Needle design facilitates one-handed operation during procedural use. Etched tip enhances visualisation during ultrasound imaging.
- Precise, rapid-firing spring mechanism provides intact core tissue sample.
- Bevelled point stylet permits easy penetration into the specimen with less trauma to surrounding tissue.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.



ORDER No.	Remarks
V-QC-20-15.0-10T	20 gauge bevelled point stylet, 15.0cm length with 10mm specimen notch.
V-QC-20-9.0-20T	20 gauge bevelled point stylet, 9.0cm length with 20mm specimen notch.
V-QC-18-15.0-10T	18 gauge bevelled point stylet, 15.0cm length with 10mm specimen notch.
V-QC-18-9.0-20T	18 gauge bevelled point stylet, 9.0cm length with 20mm specimen notch.
V-QC-16-15.0-20T	16 gauge bevelled point stylet, 15.0cm length with 20mm specimen notch.
V-QC-16-9.0-20T	16 gauge bevelled point stylet, 9.0cm length with 20mm specimen notch.
V-QC-14-20.0-20T	14 gauge bevelled point stylet, 20.0cm length with 20mm specimen notch.
V-QC-14-15.0-20T	14 gauge bevelled point stylet, 15.0cm length with 20mm specimen notch.
V-QC-14-9.0-20T	14 gauge bevelled point stylet, 9.0cm length with 20mm specimen notch.