

Endoscopy

Quantum™ Balloon Dilator



Oesophageal Balloon Dilator

DESCRIPTION

- Used to dilate strictures of the gastrointestinal tract, including strictures of the oesophagus, pylorus or colon.
- Supplied non-sterile.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.



ORDER No.	Inflated Diameter	Balloon Length	Maximum Inflation Pressure
V-QD-12x8-E	12mm/36 French	8.0cm	90psi/6.0atm
V-QD-14x8-E	14mm/42 French	8.0cm	60psi/4.0atm
V-QD-16x8-E	16mm/48 French	8.0cm	60psi/4.0atm
V-QD-18x8-E	18mm/54 French	8.0cm	50psi/3.5atm

Medflator® Inflation Device

DESCRIPTION

- Illuminated gauge for easy reading in darkened conditions.
- Fast and easy priming.
- Large threading for minimal turns.
- Ergonomic trigger allows both gradual and rapid pressure release.
- Steering handle - torque device.
- Insertion tool.
- Intended for VETERINARY USE ONLY.



ORDER No.	Remarks
V-E-T1459	Medflator® inflation device.