


Miscellaneous

Pericardial Drainage Catheter Set



DESCRIPTION

- Used to drain fluid from the pericardial sac.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.


ORDER No.

Remarks

V-PLS-911A

9 French, 20cm long catheter with six distal sideholes, an 18 gauge, 7.0 long introduction needle and a 0.038" inch diameter, 50cm long wire guide.

Martin Direct Puncture Pericardiocentesis Set



DESCRIPTION

- A direct puncture access catheter for introduction of a drainage catheter into the pericardial sac.
- The hollow introduction needle assists confirmation of correct needle access.
- A mark at 1.0cm proximal to the most proximal sidehole assists with correct catheter placement.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

ORDER No.

Remarks

V-TPT-1020U-MARTIN

10 French, 40cm long polyurethane catheter with 14cm long introduction needle.

REFERENCE

Mike Martin MVB, DVC, MRCVS, RCVS. Kenilworth, United Kingdom.