

## Ear Lavage Catheters

### DESCRIPTION

- A soft atraumatic catheter for lavage of the ear canal.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.



### ORDER No.

### Remarks

V-CELC-700U

7 French, 20cm long polyurethane catheter with four (4) sideholes.


V-CELC-850U

8.5 French, 20cm long polyurethane catheter with four (4) sideholes.

## Nasolacrimal Flushing Cannula

### DESCRIPTION

- A stainless steel cannula with matching stylet for the injection of solutions into the nasolacrimal duct.
- A stop on the cannula 1.0cm back from the tip aids in correct positioning.
- Supplied sterile in peel-open packages.
- Intended for VETERINARY USE ONLY.



### ORDER No.

### Remarks


V-SDD-23-6.0

23 gauge cannula and stylet, 6.0cm long with metal Luer lock hub.

## Nasolacrimal Flushing Catheter

### DESCRIPTION

- A PTFE catheter for the injection of solutions into the nasolacrimal duct.
- The proximal portion of the catheter is 3.0 French with the distal 1.0cm tapered to 2.0 French.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.



### ORDER No.

### Remarks

V-NLFC-3.0T-2.0

3.0 French PTFE catheter, 2.0cm long with 1.0cm distal portion, 2.0 French and plastic Luer lock hub.