



Single Line with filter vented chamber, rotating adapter and injection port



Flexible Spiral Tubing

DESCRIPTION

- Used for the infusion of fluids into the venous circulation.
- Rotating adapter is designed to prevent tangling of lines during infusion therapy.
- Supplied sterile in peel-open packages.
- Intended for ONE TIME USE.
- Intended for VETERINARY USE ONLY.

| ORDER No. | Remarks |
|-----------------------|---|
| V-FACT-14-40-FVC-1-RA | 14 French, 40cm long administration set with vented spike chamber, roller clamp and a rotating adaptor with injection port. |
| V-DPT-9.5-250-P-FM-FS | 9.5 French Flexible Spiral connecting tube set. Uncoiled length 250cm. |