

Accessories



STOPCOCKS

ORDER No. **Remarks**

V-POWSC-FLL-MLL-LP One-way plastic stopcock, female Luer lock to male Luer lock.



ORDER No. **Remarks**

V-PTWSC-2FLL-MLL-R-LP Three-way plastic stopcock, two female Luer locks to male Luer lock with right side arm.



MULTI-PURPOSE TUBING ADAPTORS

ORDER No. **Remarks**

V-PFLLA-VTA Female Luer lock to tapered catheter end. Fitting accepts 10-20 French inside diameter.



ORDER No. **Remarks**

V-PFLLA-VTA-L Female Luer lock to large tapered catheter end. Fitting accepts 20-40 French inside diameter.



ORDER No. **Remarks**

V-PMLLA-VTA-L Male Luer lock to large tapered catheter end. Fitting accepts 20-40 French inside diameter.



ORDER No. **Remarks**

V-P412 4.0-12mm double-ended tapered connector.



ORDER No. **Remarks**

V-P79 7.0-9.0mm double-ended tapered connector.

MULTI-PURPOSE TUBING ADAPTORS CONT'...

FIXATION DISC

ORDER No. Remarks

V-FD-100 Silicone fixation disc used to stabilise percutaneous drainage catheters.



INJECTION PORTS

ORDER No. Remarks

V-VIP-200 Luer lock cap to allow injection of solutions through cap.



CLAMPS

ORDER No. Remarks

V-CLAMP-W Clamp to occlude drainage catheter tubing.



ORDER No. Remarks

V-CLAMP-W-L Large clamp to occlude drainage catheter tubing.



SILICONE ADHESIVE

ORDER No. Remarks

V-ADH-2 Single tube of silicone adhesive. Recommended for repair of balloon cuffs.



LUER LOCK PLUG

ORDER No. Remarks

V-PFLLP-100 Luer lock plug to occlude female Luer lock.

