



Non-Cuffed Wire Reinforced Endotracheal Tubes – Silicone



DESCRIPTION

- The transparent silicone tube allows easy visual inspection.
- The metal spiral wire enables the position of the tube to be easily confirmed through x-ray control.
- This special design prevents the tube from kinking whilst it is in use.
- The bevel tip allows easy insertion.
- Intended for VETERINARY USE ONLY.



ORDER No.	Tube I.D.(mm)	Tube O.D.(mm)	French Size	Length (cm)
V-PATR-25	2.5	4.0	12	14
V-PATR-30	3.0	4.7	14	16
V-PATR-35	3.5	5.3	16	18
V-PATR-40	4.0	6.0	18	20
V-PATR-45	4.5	6.7	20	22
V-PATR-50	5.0	7.3	22	24

Note: Intended for use in conjunction with V-MIS malleable stiffening stylet (see below).

Malleable Stylets

ORDER No.	Remarks
V-MIS-3.0-12.5	Malleable introduction stylet. Recommended for 2.0, 2.5 and 3.0mm I.D. non-cuffed tubes.
V-MIS-5.0-35	Plastic-coated malleable introduction stylet. Recommended for 3.0 and 4.0mm I.D. cuffed endotracheal tubes.
V-MIS-12-45	Malleable introduction stylet. Recommended for 5.0, 6.0 and 7.0mm I.D. cuffed endotracheal tubes.