



## Cuffed Wire Reinforced Endotracheal Tubes – Silicone



### DESCRIPTION

- The transparent silicone tube allows easy visual inspection.
- The metal spiral wire enables the position of the tube to be easily confirmed through x-ray control.
- This special design prevents the tube from kinking whilst it is in use.
- Low pressure cuff enables atraumatic positioning.
- The bevel tip allows easy insertion.
- Intended for VETERINARY USE ONLY.




ORDER No.	Tube I.D.(mm)	Tube O.D.(mm)	French Size	Length (cm)
V-PETR-6	6.0	8.7	26	28
V-PETR-7	7.0	10.0	30	30
V-PETR-8	8.0	11.3	34	32
V-PETR-9	9.0	12.7	38	34
V-PETR-10	10.0	14.0	42	36

**Note:** Intended for use in conjunction with V-MIS malleable stiffening stylet (see bottom of page 5).

## Minitracheostomy Tube

### DESCRIPTION

- The polyurethane minitracheostomy tube is supplied with a blunt obturator and airway adaptor.
- The blunt obturator allows atraumatic placement.
- Removable airway adaptor allows the introduction of suction catheters into the bronchial airways for aspiration.
- A small cap is fitted to allow the tube to be sealed.
- Intended for VETERINARY USE ONLY.
- Intended for ONE TIME USE.



ORDER No.	Remarks
V-CMST-400-BO	4.0mm I.D. minitracheostomy tube with obturator. Intended for one-time use.