


Endotracheal Tubes with Replacement Accessories



DESCRIPTION

- A cuffed all silicone endotracheal tube used to deliver anaesthesia gases after positioning by oral intubation into the desired position.
- The balloon cuff can be inflated to prevent gas loss during assisted respiration.
- The proximal fitting of tubes up to 12mm accepts a 15mm adaptor. A funnel type silicone fitting is supplied with 14 and 16mm tubes.
- Repair balloon cuffs and inflation lines are available as well as tubes of silicone adhesive.
- Intended for VETERINARY USE ONLY.



REPLACEMENT ACCESSORIES

ORDER No.	Tube I.D. (mm)	Tube O.D. (mm)	French Size	Length (cm)	Balloon Length (cm)	Inflation Lines	Repair Cuff	Airway Adaptor	Silicone Adhesive
V-PET-2	2.0	3.7	10	16	2.0	V-RIL-S	V-PRK-2	V-P309-2	V-ADH-2
V-PET-3	3.0	4.7	14	16	2.0	V-RIL-S	V-PRK-3	V-P309-3	V-ADH-2
V-PET-4	4.0	6.0	18	18	2.0	V-RIL-S	V-PRK-4	V-P309-4	V-ADH-2
V-PET-5	5.0	7.3	22	20	3.0	V-RIL-S	V-PRK-5	V-P309-5	V-ADH-2
V-PET-6	6.0	8.7	26	28	4.0	V-RIL-S	V-PRK-6	V-P309-6	V-ADH-2
V-PET-7	7.0	10	30	30	5.0	V-RIL-S	V-PRK-7	V-P309-7	V-ADH-2
V-PET-8	8.0	11	33	32	5.0	V-RIL-S	V-PRK-8	V-P309-8	V-ADH-2
V-PET-9	9.0	12.3	37	34	5.0	V-RIL-S	V-PRK-9	V-P309-9	V-ADH-2
V-PET-10	10	14.3	43	36	5.0	V-RIL-S	V-PRK-10	V-PRK-10	V-ADH-2
V-PET-11	11	15.7	47	40	5.0	V-RIL-S	V-PRK-11	V-P309-11	V-ADH-2
V-PET-12	12	17	51	40	5.0	V-RIL-S	V-PRK-12	V-P309-12	V-ADH-2
V-PET-14	14	19	57	48	7.5	V-RIL-S	V-PRK-14	V-P309-14	V-ADH-2
V-PET-16	16	22	66	70	8.0	V-RIL-L	V-PRK-16	V-P309-16	V-ADH-2

