



DESCRIPTION

- A closed-end, single lumen silicone catheter used for the administration of feeding solutions.
- Catheter can be introduced orally or nasally into the stomach or intestine.
- Lightweight low profile proximal fitting.
- Intended for VETERINARY USE ONLY.

ORDER No.	Remarks
V-PFT-3.5-20	3.5 French nasogastric Feeding tube, 20cm long. Closed endhole design.
V-PFT-5-20	5 French nasogastric feeding tube, 20cm long. Closed endhole design.
V-PFT-5-35	5 French nasogastric feeding tube, 35cm long. Closed endhole design.
V-PFT-6-35	6 French nasogastric feeding tube, 35cm long. Closed endhole design.
V-PFT-8-55	8 French nasogastric feeding tube, 55cm long. Closed endhole design.
TV-PFT-10-55	10 French nasogastric feeding tube, 55cm long. Closed endhole design.

Wire Guide Stiffening Stylets

ORDER No.	Remarks
V-SF-28-60	Wire guide stylet for V-PFT-6-55.
V-SF-38-60	Wire guide stylet for V-PFT-8-40 and V-PFT-8-55.
V-SF-45-60	Wire guide stylet for V-PFT-10-40 and V-PFT-10-55.