



AQUACEL
Hydrofiber® Wound Dressing



ben^{son} medical industries inc.

Research Supplies and Equipment



Benson Medical Industries Inc. is a Canadian owned and managed company, specializing in the sale and service of veterinary, hospital and dental products since 1974.

Benson Medical's head office is located in Markham, Ontario. We also serve our customers from branch offices and sub-distributors located across Canada, as well as designated sales representatives, office support staff and Clinical Product Specialists.

Benson Medical provides technical support and maintains a fully equipped service facility with factory trained personnel for both in-shop and field service repair.

Benson Medical Industries carries many items to meet your research needs. If you do not find what you are looking for, please contact us for assistance.

Why buy “Canadian”?



Benson Medical has adjusted and will continue to adjust our pricing to mirror the changes in the US/Canadian exchange rate. However, there are a number of other important factors that the purchaser should consider when comparing a product purchased from Benson Medical versus an American counterpart.

You may want to ask the following questions:

Are there any differences in the product warranty or support?

***We fully support what we sell and offer in-service training for our customers.**

Can we obtain rental product in the event of an equipment problem?

***We have some rentals available to locations that purchase their equipment from us.**

Does the equipment pass Canadian CSA electrical standards?

***We only sell CSA approved equipment.**

What exchange rate do I get if I order directly from the US?

***All of our prices are in Canadian funds.**

What are my shipping, brokerage and/or insurance costs when ordering equipment from the US?

***Like the CSA electrical approvals, these importation costs are included in our prices.**

If our Canadian list prices are slightly higher than those quoted from a US-based company, we hope you consider all of these associated costs and benefits before making your decision. Benson Medical offers a personalized service that is geared specifically for the Canadian market. We value your business.

Contents

Anaesthesia Machines and Accessories.....	4-9
Tabletop	4
Wall Mount	5
Mobile/Wheeled Base.....	6-7
Multi-Patient/Pre-Set	8
Accessories.....	9
Vaporizers & Fillers	9
Waste Gas Management.....	10
Ventilators	11
Breathing Circuits	12
Face Masks.....	13
Lighting.....	14
Endotracheal Tubes, Laryngoscopes	14
V-Gel®	15
Rodent Intubation & Ventilation Supplies.....	16
Adaptors, Resus. Bags & Rebreathing Bags	17
Patient Warming & Surgical Instruments	18-19
Induction Chambers	20
Oxygen Supplies & Accessories	21
Scales	22
SpO ₂ and EtCO ₂ Monitors.....	22-23
Blood Pressure Monitors	24
Mindray Monitors.....	25
Midmark & Cardell Monitors	26
Bionet Monitors	27
Microscopes & IV Pumps	28
Autoclaves & Ultrasonic Cleaners.....	29
Stainless Steel.....	30
Imaging Accessories & Euthanasia	31
Consumables	32-33
Aquacel	33
Prevail/Ogena Cleaning Products.....	34-35
Find Your Local Rep.....	35

Anaesthetic Machines and Accessories

Table Top—Non Rebreathing



Vaporizer Not Included

'PAM' Portable
Anaesthesia Machine

[#07-8914722](#)

MRI compatibility (Option)



Vaporizer Not Included

M3000
Anaesthesia Machine

[#M3000](#)

MRI compatibility (Option)



V-1 Table Top
Anaesthesia Machine

[#901801](#)



Vaporizer Not Included

VMR Table Top Assembly

[#94805059](#)



Vaporizer Not Included

Tabletop Research Anaesthesia Machine:
200cc-4 LPM O₂ Flowmeter

[#20050](#)

Other configurations available



TT2: Tabletop
Anaesthesia Unit

[#TT-17000](#)

Table Top—Rebreathing



Vaporizer Not Included

M4000 Tabletop
Anaesthesia System

[#M4000](#)



Vaporizer Not Included

VME Tabletop
Anaesthetic Machine

[#91800121](#)



Vaporizer Not Included

CDS2000 Tabletop
Anaesthetic Machine

[#20027](#)

Wall Mount—Non Rebreathing



V9 Wall Mount

[#901802](#)



Vaporizer Not Included
'PAM' Wall Mount

[#07-891472*](#)

*Specify wall mount



Vaporizer Not Included
Basket style may vary

M2350 Wall Mount

[#M2350](#)

Wall Mount—Rebreathing



V9A Wall Mount

[#901803](#)



Vaporizer Not Included

VMS Wall Mount

[#91800166](#)



Vaporizer Not Included

M2300 Wall Mount

[#M2300](#)



Vaporizer Not Included

M2500 Wall Mount

With Articulating Arm

[#M2500](#)



Vaporizer Not Included
Versa II Wall Mount

[#078910545](#)

Wall Mount Bracket

[#078914013](#)



Vaporizer Not Included

CDS2000 Tabletop

[#20027](#)

Wall Mount Kit

[#22029](#)

Mobile—Non Rebreathing



V10 - Mobile Anaesthetic Machine

#901804

Options:

Single & Double E-tank Mount, for O₂ & CO₂

Dual Vaporizers

Carrier Gas Options



Compact Convertible Anaesthesia System

#07-8915662

Options:

E-Tank Mount

Mayo Tray

Vaporizer Not Included



Vaporstick with 0.2 to 4LPM Flowmeter

#20036

0-1000mL Flowmeter

#20037

Options:

Dual O₂ E-tank Mount

Mayo Tray



M1050 Anaesthesia Machine

#M1050

Options:

Carrier Gas Options

E-tank Mount (shown, not included)

Oxygen Concentrator

Vaporizer Not Included
E-Tank Not Included
Basket style may vary



M1000 Anaesthesia Machine

#M1000

Options (some shown, not included):

Dual O₂ E-tank Mount

Stainless Steel Tray

Carrier Gas Options

VMR Mobile Anaesthetic Machine

With 0-1000mL flowmeter

#91800244

With 0-4LPM Flowmeter

#91800250

Options:

E-Tank Mount



Vaporizer Not Included

Vaporizer Not Included
E-Tanks Not Included
Basket style may vary



M6100 Anaesthesia Machine with Oxygen Concentrator

#M6100

Options:

IV Pole

Carrier Gas Options

Versa II, Non-Rebreathing Anaesthesia System

#07-8915660

Options:

Dual Vaporizers



Vaporizer & E-tank - Not Included

Mobile—Rebreathing



V10A - Mobile Anaesthetic Machine with CO₂ Absorber

#901805

Options:

Single & Double E-tank Mount, for O₂ & CO₂

Dual Vaporizers

Carrier Gas Options

Bain Arm



M1200 Anaesthetic Machine

#M1200

Options (some shown, not Included):

E-Tank Mount

Carrier Gas Options

Bain Arm

Vaporizer Not Included
E-Tanks Not Included



VMS Standard Anaesthetic Machine

#91800070

Options:

E-Tank Mount

Shelf



CDS9000

#20021

Options:

Dual flowmeter

Shelf & Drawer

E-Tank Mount



M2000 Anaesthetic Machine

#M2000

Options (some shown, not included):

E-Tank Mount

Carrier Gas Options

Dual Vaporizer

Circuit and Bag Organizer

Bain Arm

Vaporizer Not Included
Oxygen Tanks Not Included



VME 2 Mobile

#91805112

Options:

E-Tank Mount

Dual Flowmeter

Dual Vaporizer

Shelf



Versa II
Anaesthesia System

#07-8910545

Options:

Dual Vaporizers



M6000 Anaesthesia Machine with Oxygen Concentrator

#M6000

Options:

Bain Arm

Carrier Gas Options

Vaporizer Not Included
E-Tank Not Included

Multi-Patient & Pre-Set



RC² - Rodent
Circuit Controller

[#922100](#)

Mobile Unit

[#922120](#)



RC²⁺ - Rodent
Circuit Controller Plus

[#905300](#)



SAS-3 - Smart
Anaesthesia System

[#07-8909689](#)



IMPAC⁶ - Integrated Multi-Patient
Anaesthesia Center

[#901808](#)



Link7 Anaesthesia &
Evacuation System

[#07-8915712](#)



Compac⁵ Mobile Unit

[#901812](#)

[#901816](#)




Multi-Station Lab
Animal Machine:

2-Station: [#20300](#)

4-Station: [#20302](#)

6-Station: [#20304](#)



Multiple Rodent Surgical Board

[#RES3200](#)



AMD-3

[#AMD-17004](#)



Dual Diverter Manifold

[#07-8914724](#)

Tri-Diverter Manifold

[#07-8914725](#)



Dual Anaesthesia
Set:

[#20382](#)



Multiplex Delivery System:

2 Station: [#07-8914732](#)

3 Station: [#07-8914733](#)

4 Station: [#07-8918185](#)



Balanced Flow
Multi-Station Manifolds:

2 Port: [#RES535](#)

4 Port: [#RES536](#)

6 Port: [#RES537](#)

Accessories



Universal Control Arm
[#07-8914012](#)



Mask Stabilizer
[#07-8909415](#)



Bain Circuit Adaptor
[#MA2013](#)



Humidifier for Rodent NRB System
[#07-8914735](#)



High Pressure Alarm
[#V7320](#)



Dual E-tank Manifold
 Available for Many Units



Circuit Secure
[#RES520](#)



Mayo Tray/Shelf
 Available for Many Units

Vaporizers



New TEC III Isoflurane, Cagemount, Key-Fill
[#53-T3ISO-NVKL](#)

New TEC III Isoflurane, Cagemount, Funnel-Fill
[#53-T3ISO-FF](#)



Sevoflurane options available

All vaporizers listed above are available with a selecatec mount as an option.



Vapofil Isoflurane key-filler
[#12016](#)



Isoflurane anti spill adaptor
[#51030](#)



Waste Gas Monitoring Badge and Analysis Report
[#931443](#)

Waste Gas Scavenging



AES Active Evacuation System

#20500



EVAC 2 Scavenging System

#07-8918181



Extract-All Bench Evac Unit

#933101



Extract-All Cabinet Evac Unit

#933105



Filter Mate Portable Extractor

#935101



Lab Animal Machine Evacuation System

#20504



SI-6 Scavenger Interface Device

#931445



Scavenging Cube

#931600



F/Air Canister—8 pack

#20522



VaporGuard Charcoal Filter 8 pack

#931401



Charcoal Pre-Pack

#20510



19mm Scavenge Tubing

#10024



Scavenger Interface Valve

#931439



Gas Evacuation for Central Vacuum

#EVC629



VaporGuard XL-R Refillable Charcoal Filter (2 pack)

#931402



Bulk Carbon Container

#101500



Scoop Stereotaxic Waste Gas Remover

#078914736



SOMNI Catch

#WG-15005



Gas Evacuation for Fan or Blower

#EVC630



WAG Canister - 4 pack

#07-8909457

Ventilators



Anaesthesia Work Station (AWS)
Small Animal Ventilator

#70008

200mL Attachment

#70010



Small Animal Ventilator

#SAV04-SYSTEM



Merlin Small Animal Ventilator

#VSL10-1002-02



SAV 2500 Ventilator

#V725000

MRI 2250 Ventilator

#V727000



Model 2002Pro Ventilator

#000A5766

Adjustable I:E Ratio

#000A5770

MRI compatibility (Option)



Model 3002 Ventilator

300-1600mL

#91806001



Rat Ventilation
Nosecone

#921465



ADS2000 Anaesthesia Delivery
System w/12 hour battery

#2000-2000

Stand with Vaporizer Clamp

#1000-5000



Nuffield 200 Ventilator

#5007286



Ventilator bellows options:

0-300mL

300-1600mL

1600-3000mL


Breathing Circuits



60" Bain Circuit
[#BD2491](#)



BAIN (Blue Inner Tube)
 60" Length
[#50521](#)
 40" Length
[#50520](#)



Ultra Bain System, Deluxe with Delrin & Scavenger Adaptors
[#50515](#)
 Ultra Bain Circuit with Scavenger Adaptors
[#50516](#)
 (Re-breathing Bag not Included with Circuit)



UniFlow Nosecone
[#EQFU001](#)



Mapleson D - Rodent NRB
[#07-8909636](#)



NRB Circuit with Nosecone
 9mm nosecone: [#921461](#)
 12.5mm nosecone: [#921462](#)
 14mm nosecone: [#921463](#)



Mouse and Rat NRB
[#RES500](#)



Modified Jackson Rees
 (without bag)
[#50510](#)



60" Universal 'F'
[#50523](#)
 40" Universal 'F'
[#50522](#)




Fixed 'Y' Piece
 60" Length
[#50502](#)
 40" Length
[#50500](#)

Pediatric 'Y', Fixed, 60"
[#50508](#)
 Pediatric Circle Circuit (MRI)
 25 foot length
[#50508-25](#)
 30 Foot length
[#50508-30](#)

Custom lengths available

Animal Face Masks

Autoclavable Animal Masks



Animal Mask, Sz 0

[#50128](#)

Animal Mask, Sz 2

[#50130](#)

Animal Mask, Sz 4

[#50132](#)

Animal Mask Set

(One of each Size)

[#50134](#)

Animal Mask, Sz 1

[#50129](#)

Animal Mask, Sz 3

[#50131](#)

Animal Mask, Sz 5

[#50133](#)

Rodent Face Masks



Small Mask (up to 150g)

[#305](#)

Large Mask (200-350g)

[#307](#)

Medium Mask (150-200g)

[#306](#)

Extra Large Mask (350g+)

[#308](#)

Replacement O-rings and 12 mm Latex Sheets available

Nose Cone, 9mm, Mouse

[#921609](#)

Nose Cone, 12mm,

Large Mouse to Medium Rat

[#921612](#)

Nose Cone, 14mm, Rat

[#921614](#)

Rodent Mask (no diaphragm)

[#20350](#)

Rodent Mask Diaphragm

[#20352](#)

Stereotaxic Mask:

Mouse

[#921435](#)

Rat

[#921434](#)

Rat-Taxic:

[#078917519](#)

Finger Cot Nosecone Kit

[#921610](#)

Replacement Finger Cot
Black Latex

[#190501](#)

Lighting



LED 120, 130



255



Celestial Star

LED & Halogen Surgical & Exam Lights

Ceiling, Wall and Mobile Options

MRI Versions Available



E.T Tubes, Laryngoscopes & V-Gel®

PVC Endotracheal Tubes:

Sizes 2.0-2.5mm, uncuffed

[#50800-size](#)

Sizes 3.0-10.0mm in 0.5mm increments

[#50802-size](#)

Set of 15 (3.0-10.0mm)

Size 11

Size 12

Size 13

[#50250](#)

[#50252](#)

[#50254](#)

PressureEasy ET
Cuff Pressure
Controller

[#21800](#)



LED Laryngoscope Handles:

F/L Green LED Handle, Stubby "AA"

[#19060LED](#)

F/L Green LED Handle, Penlite, "AA"

[#19062LED](#)

F/L Green LED Handle, Medium "C"

[#19064LED](#)



Silicone Endotracheal Tubes:



Sizes 2-14mm in 1mm increments

[#ET\(size\)](#)



Crystal



Macintosh



Miller

Penlon Laryngoscopes:

Premier, Diamond and ISO Red laryngoscopes with standard, stubby or penlight handles and Macintosh, Miller, Seward, Robertshaw, McCoy or Wisconsin blades.

MRI compatible and disposable options available.

Ask your rep for order numbers and prices.



Seward



Robertshaw



RABBIT v-gel® Sizing Guide

Size	Ideal body weight of rabbit	Product code
R1	0.6 – 1.5kg (1.3 – 3.5lb)	D30001
R2	1 – 2kg (2 – 4.5lb)	D30002
R3	1.8 – 3.5kg (4 – 8lb)	D30003
R4	2.5 – 4kg (5 – 9lb)	D30004
R5	3.5 – 5kg (8 – 11lb)	D30005
R6	4.5kg+ (10lb+)	D30006



CAT v-gel® Sizing Guide

Size	Ideal body weight of cat	Product code
C0	0.5 – 1kg (1 – 2lb)	D40000
C1	1 – 2kg (2.5 – 4.5lb)	D40001
C2	1.5 – 3.5kg (3 – 8lb)	D40002
C3	3 – 5kg (6.5 – 11lb)	D40003
C4	4.5 – 6kg (10 – 13lb)	D40004
C5	6kg+ (13lb+)	D40005

V-Gel® Advanced Features:

- Fast, easy, safe and stress free insertions
- Anatomical matching features combined with a soft material to give a high quality pressure seal
- No post-operative coughing or gagging
- Low airway breathing resistance due to the large airway channel within the device
- Suited to ventilated as well as spontaneously-breathing patients
- Greater seal pressure – even better for ventilating
- High quality pressure seal restricting leakage of volatile anaesthetic agents, thus improving health and safety for staff
- Tip inserts into the upper esophagus to help prevent and manage regurgitation of stomach contents
- Bowl now envelops the whole larynx and epiglottis sits inside – avoids “blue tongue” in rabbits!
- More space and security for dental procedures
- Safer, smoother recoveries

<http://docsinnovent.com/>



D-grip® tube holder is a fully adaptable airway management support system for supporting the load borne by anaesthetic breathing circuits, airway devices, facemasks, monitoring cables and sampling lines. It is the perfect complement to v-gel® and offers a wide range of positioning options to optimise airway management.

#D10013



VetLube
30mL Spray Bottle

#D14002



Rodent Intubation & Ventilation Nosecones



Rodent Tilting
WorkStand

#70018



Mouse Intubation Pack

#71002

(Otoscope and Lidocaine Jelly
not included)



Short Mouse Specula,

Autoclavable

#71030S



Rat Intubation Pack

#71000

(Otoscope and Lidocaine Jelly
not included)



Long Mouse Specula,

Autoclavable

#71030



Mouse Lidocaine Applicator

#71038



Rat Lidocaine Applicator

#71060



Mouse ET Tube Introducer

#71036



Rat Specula,

Autoclavable

#71050



Rat ET Tube Introducer

#71058



Mouse ET Tube, 1.22mm

5 pack

#71034



Rat Ventilation Nosecone

#921465



Rat ET Tube, 1.27mm

5 pack

#71032

Adaptors and Connectors

6mm Male x 15mm Male tapered fresh gas outlet adaptor

#50305



23mm MJR tapered vaporizer outlet adaptor

#50382



22mm Male/15mm Female x 22mm Male/19mm Female

#50316



19mm Male x 22mm Female

#50304



Vaporizer Connector Outlet Adaptor

#50632




Vaporizer Connector Inlet Adaptor

#50622

19mm Male x 30mm Female (for some styles of pop off valves)

#KAB-1930D



22mm Male x 22mm Female Soft plastic sleeve

#836-113-100F



Scavenger "Y" Adaptor (19 mm M x 19 mm F x 19 mm F)

#50376



Circuit Elbow Without Port

#255-3507-800

Circuit Elbow With Port

#255-3523-804



NON re-breathing Adaptor (For MJR circuit)

#20700V



Wye Piece Adaptor (7mm ends)

#30614-55



DELIRIN ADAPTOR

(22 mm Male - 15 mm Female x 23 mm Male - 15 mm Female)

#50518



Male Adaptor - Quick Connect

#850-0030-000



CP – Female Quick Connect

#850-PLC17004



Induction Chamber Fresh Gas Connector

#960-0030-000



Re-Breathing & Resuscitation Bags



Latex Free

Re-Breathing Bags:

500cc

3 litre

Easy Grip Adult, Autoclavable Silicone

#10520B

#10530B

#01BM3600

1 litre

4 litre

Easy Grip Child, Autoclavable Silicone

#10524B

#10532B

#01BM3610

2 litre

5 litre

Easy Grip Infant, Autoclavable Silicone

#10528B

#10534B

#01BM3620



Easy Grip BVM Resuscitator, Large

#4000-BL (pictured)

Easy Grip BVM Resuscitator, Medium

#4010-BL

Easy Grip BVM Resuscitator, Small

#4020-BL

Patient Warming



Heat Therapy (water) Pump

#70308

*Extension tubing available *



Fluid Warmer with
LCD Display

#FW-300VET




Heated Rodent Surgical Table
For Active Scavange

#RES4000




Heated Rodent Surgical Table
For Passive Scavange

#RES4100



DSx Vented Warming Table

#965500



ProStation - Vented Warming Table

#07-8914731



Blue Vinyl Pad 12" x 18"

#VRP1218

Blue Vinyl Pad 22" x 30"

#VRP2230

Blue Vinyl Pad 24" x 60"

#VRP2460



Fleece Cover for 14" x 22" Pad

#FC1218

Fleece Cover for 22" x 30" Pad

#FC2230

Fleece Cover for 24" x 60" Pad

#FC2460



Small Black Vinyl Pad 16" x 20"

#V016

Medium Black Vinyl Pad 20" x 29"

#V029

Large Black Vinyl Pad 20" x 44"

#V044



Disposable Warming Pads, case of 20

Neck Collar - 4" x 19"

[#ST-003](#)

Small - 13.5" x 18.5"

[#ST-017](#)

Medium - 16" x 22"

[#ST-020](#)

Large - 20" x 36"

[#ST-024](#)



Heated Hard Pad 6" x 8"

[#70305](#)

Heated Hard Pad 9" x 12"

[#70304](#)

Replacement Connector Kit

[#000A4083](#)

Heated Hard Pad 15" x 24"

[#70302](#)

Heated Hard Pad 21" x 24"

[#70300](#)



Hot Dog® Controller

[#WC71V](#)

Hot Dog® Multi-Port Controller

[#WC77V](#)

Extra Small Blanket 10" x 12"

[#V101](#)

Small Blanket 10" x 23"

[#V102](#)

Medium Blanket 16" x 22"

[#V103](#)

Large Blanket 22" x 31"

[#V104](#)

X-Large Blanket 27" x 47"

[#V106](#)

Electrical Cable, 4m

[#A101](#)

Wire Storage Rack for Blankets

[#A108VET](#)



Instruments



High Quality Stainless Steel Instruments:

Dental

Surgical

Left Handed

Micro Instruments

Call for more information

Induction Chambers



Induction Chamber, 1L - 3.75"w x 4.5"d x 3.75"h

[#941443](#)

Induction Chamber, 2L - 3.75"w x 9.0"d x 3.74"h

[#941444](#)

Induction Chamber, 7L - 8.25"w x 8.75"d x 5.75"h

[#941448](#)

Induction Chamber, 15L - 15"w x 7.5"d x 7.7"h

[#941454](#)

Replacement gaskets available



Induction Chamber 1.5L - 7.5" x 5.5" x 4"

[#93805106](#)

Induction Chamber 5.25L - 11" x 5.75" x 7"

[#93805107](#)

Induction Chamber 18.5L - 16" x 8" x 10.5"

[#93800038](#)

Induction Chamber 45L - 20" x 12" x 12"

[#20220601](#)



Sliding Top Chamber:

X-Small (1.25 L):

[#078917855](#)

Small (3.0 L):

[#078933388](#)

Large (5.4L):

[#078933387](#)



Red Sliding Top Chamber:

Small - 4.25" x 4.25" x 4.25"

[#IR-3005](#)

Medium - 9.25" x 4.25" x 4.65"

[#IR-3006](#)

Large - 11.75" x 6" x 5.75"

[#IR-3007](#)



Mini Chamber 10" x 5.5" x 5.5"

[#55205](#)

Small Chamber 16" x 7" x 7"

[#55212](#)

Larger sizes available



Induction Chamber, 1 L, vented

[#942103](#)

Induction Chamber, 2 L, vented

[#942102](#)

Induction Chamber, 7 L, vented

[#941447](#)



Vapor-Vac Hooded Chamber (9.0 L)

[#78933385](#)



Active Induction Chamber

[#IA-3009](#)



Red Vented Chamber

Medium - 9.25"x 4.25" x 5.75"

[#IR-3010](#)

Large - 11.75" x 6" x 5.75"

[#IR-3011](#)

Flow Meters, Concentrators, Regulators & Accessories



Oxygen Regulator - 'E' yoke
With Male DISS Fitting

#50601

* Also available for N₂O & Nitrogen*



Oxygen Regulator for Large Tanks
With Male DISS Fitting

#50600

* Also available for CO₂, N₂O & Nitrogen*



0-25 LPM, Barb Outlet, with Gauge

#01RO01003-A



Oxygen Concentrator

#OC-21002



Oxygen tank not included

Single Cylinder Cart (D,E)

#6110

Dual Cylinder Cart (D,E)

#WE-HCC5-M

12 Cylinder Cart (D,E tanks)

#MCC1



0-15 LPM Oxygen Flow Meter
With Female DISS - Hand Tight

#15408-03

0-12 LPM CO₂ Flow meter

#FM-803C

0-15 LPM Air Flow Meter

#15508-13



Oxygen Concentrator

#OC4000



Oxygen Cylinder Wrench

#20804



Disposable Manometer

#55356



Oxygen Adapter

#1244-3DW



Dual Oxygen Outlet (DISS)

#50932



Custom length hosing is available with various fittings (DISS, Puritan, etc.) for Oxygen, CO₂, N₂O, Nitrogen, Vacuum and Air.

Call for details.



Scales



LCD Cat Scale - Up to 44lbs
20" x 14"

[#110-5090-07](#)



Indicator floor stand sold separately

Veterinary Model PS500 Weigh Scale - Up to 550lbs

22" x 22"

[#816965001903](#)

22" x 42"

[#816965001453](#)

22" x 36"

[#816965001897](#)



Indicator floor stand sold separately

Large Animal Scale - Up to 3000lbs
39.37" W x 78.75" L x 3.93" H

[#816965004379](#)



Veterinary Model MS20S Weigh Scale - Up to 44lbs
20.5" x 11.5" x 2.25"

[#816965001002](#)

SpO₂ & EtCO₂ Monitors



Masimo RAD57 SpO₂ Monitor

[#9216](#)

Lingual probe

[#1895](#)

Transflectance probe

[#1896](#)



Masimo Emma Immediate Capnograph

[#3639P](#)



Viamed hand held Capnographs

Mainstream CO₂ / SpO₂

[#4410530](#)

Mainstream CO₂ / SpO₂ with Multigas

[#4410560](#)



Capnostream 35 EtCO₂ & SpO₂ Monitor

[#PM35-VET-KIT](#)

Small Filterline Set

[#SV-02-0338](#)

Large Filterline Set

[#01-02-0392](#)



Vetcorder Pro - SpO₂, ECG, Temp (Optional CO₂)

[#48561](#)



Vetcorder Pro Airmate CO₂ Module

[#48520](#)



2500A Handheld Pulse –Ox

#2500A Vet



Charging Stand for 2500A

#2500C-UNIV



9847V Handheld SpO₂ & CO₂ Monitor

#9847V



CO₂ Sensor for 9847V

#9840SA



Airway Adaptor Tube pkg of 3 (9847V)

9840AAT-V



LifeSense– Pulse-Ox, EtCO₂

#5775-002



RespSense– Capnograph Only

#5776-002



7500FO Pulse-Ox

#7500FO

*MRI compatible cable *

Max heart rate 321bpm



Moisture Trap & Filter

#SENS TF-1-10

Filters - pkg/25

#8517-001



CO₂ Sample Line (pkg/5)

#5821-000



Fiber Optic Cable for MRI

20ft - **#8000FC-20**

30ft - **#8000FC-30**



Fiber Optic Cable for MRI

20ft - **#8000FI-20**

30ft - **#8000FI-30**



Lingual Clip

#2000SL




Flexible Sensor

#2000SA



Transflectance Sensor

#2000T



PS900 Cable Oximeter (requires Windows based laptop/tablet)

Rat:

#78938072

Mouse:

#78938073

Veterinary :

#78938074

Blood Pressure Monitors

The new petMAP+ II offers continuous monitoring, SD card data transfer, rechargeable battery and optional SpO₂, ECG, EtCO₂ & temperature. Specifically designed for cats and dogs.

petMap+ II™ Graphic Blood Pressure Monitor, 5" Display

#R7450

petMap+ II™ Graphic Blood Pressure Monitor with Nonin SpO₂, 5" Display

#R7454

petMap G3 Graphic Blood Pressure Monitor with Temperature Option

#BM7350



ECG Module

#9026



EtCO₂ Module

#9027



34" Temperature Cable

#8040



Charger

#9043



Replacement Cuffs Available in Sizes 2 - 13



Cover for G3

#9057

Cover for +II, 5" display

#9038-1BLK



Pole Mount

#9042




Cardell Insight Diagnostic Monitor:

Cardell Insight - NIBP only

#8014-001

Cardell Insight - NIBP & SpO₂

#8015-001



MemoDiagnostic HDO Pro

#MDPRO



Parks Doppler 811B Complete

#811B



Vmed Vet Dop2 - Complete Set

#BF2S

Multi-parameter Monitors

mindray



iMEC8-Vet - 8.4" LED Display, ECG, Resp, NIBP, SpO₂, ST and Arrhythmia, IBP, Temperature capability & recorder.

[#6307F-PA00002](#)

With Mainstream CO₂:

[#6307F-PA00013](#)

With Sidestream CO₂:

[#6307F-PA00015](#)



iPM12-Vet - 12.1" LED Display, ECG, Resp, NIBP, SpO₂, ST and Arrhythmia, IBP, Temperature capability & recorder.

[#9223F-PA00004](#)

Add Mainstream CO₂:

[#115-035565-00](#)

Add Sidestream CO₂:

[#115-030779-00](#)

Add Multigas, Invasive Pressures, Cardiac Output:

[#115-030784-00](#)



SpO₂ Probe

[#9101-10-58134](#)



SpO₂ Cable

[#9101-10-58134](#)



ECG Trunk Cable

[#0010-30-42719](#)



ECG 3 Lead Wire

[#0010-30-42719](#)



ECG Clips (5 pack)

[#9101-20-58104](#)



Rolling Stand

[#045-003428-00](#)



Temperature Probe

[#0011-30-37392](#)



Thermal Paper

[#0011-30-37392](#)



Sampling Line

[#115-043017-00](#)



Water Trap

[#125-000005-00](#)



Cardell Insight - NIBP only

[#8014-001](#)

Cardell Insight - NIBP & SpO₂

[#8015-001](#)



Midmark Multiparameter Monitor with Cardell BP, SpO₂, ECG, Resp, Temp, 12.1" Touch screen & USB port. (CO₂, Printer, IBP, Multigas options available)

Multiparameter Monitor with Printer

[#8019-001](#)

Multiparameter Monitor without Printer

[#8019-002](#)

Multiparameter Monitor with IBP & Printer

[#8019-003](#)



Max-Fast Sensor

[#Maxfast-1](#)



ECG 3 Lead Wire Set

[#016-1604-00](#)



NIBP Cuffs
Sizes 1-5, 8 & 10

[#SV-'size'](#)



Flat ECG Clip
(pkg of 5)

[#CD-0019](#)



Printer Paper
4 rolls

[#PAPER-4F](#)



SpO₂ Sensor

[#V-SAT](#)

Extension Cable

[#01-02-0183](#)



Rolling Stand

[#8008-005](#)



ECG Trunk Cable
(Touch, 9500, Max-12)

[#016-1743-00](#)

For 8019:

[#016-10177-00](#)



Temperature Probe

[#590004](#)

For 8019:

[#016-10164-00](#)



IBP Cable
(Touch, Max-12)

[#650-208](#)

For 8019:

[#016-10180-00](#)



Esophageal Probe
(Touch, 9500, Max-12)

[#EP-XS, #EP-S, #EP-L](#)

For 8019:

XS - [#016-10232-00](#)

SM - [#016-10233-00](#)

LG - [#016-10234-00](#)



IBP Transducer

[#016-1586-00](#)

With Admin Kit

[#016-1587-00](#)



Mainstream CO₂ Sensor Module

[#CSTAT-5](#)



Sidestream CO₂ Sensor Module

[#LoFlo](#)



We carry all Cardell accessories, not all are listed



BM3 Pro - SpO₂, NIBP, RR, Temp, ECG, BT-Link & BT-link mobile software, printer, Wifi, USB memory & 8" touch screen

#BM3VetPro

BM5 Pro- As above with 2 - Temp, optional 2 - IBP (10.4" touch screen)

#BM5VetPro

BM7Pro - As above with 2 - IBP, 2 - USB, SD card slot (12.1" touch screen)

#BM7VetPro

EtCO₂ or Dual gas can be added to any of the above monitors



Eosophageal Probe, ECG & Temp:
L:45mm, D:6mm: **#B-EPRB-XS**
L:70mm, D:8mm: **#B-EPRB-M**



Transflectance Sensor
#B-SPO2-R



SpO₂ Y Sensor
#B-SPO2-V



SpO₂ Extension Cable
#B-SPCBL-N



3 Lead ECG Wire Set
#B-WIRE-N



3 Lead ECG Trunk Cable
#B-CBL-N



Temperature probe
#TEM-PRV-R



Regular Airway Adaptor
#6063-00



Large Airway Adaptor
#3473ADU-00



Small Airway Adaptor
#3473INF-00



Dual Gas Module
#DG-SENSOR

Rolling Stand for BM3, 5 & 7

#BM-CART-H

Microscopes, Scopes & Video Scopes



Model 215 - Has a Seidentopf head for comfortable viewing, coaxial focusing, low drive mechanical stage and LED illumination with conveniently located intensity control. LED light rated up to 50,000 hours.

#215-RLED-ASC

The M17 compound microscope series features infinity corrected optics, ideal for college laboratories, veterinary and medical applications. Supplied with binocular head, this research level microscope also incorporates a halogen illumination system and wide field 20mm eyepieces for a very crisp image.

Infinity Corrected Binocular Microscope:

#M17B-P

Infinity Corrected Trinocular Microscope:

#M17T-P



New additions to the M10 line feature a high-resolution 3 megapixel digital camera for capturing live images of up to 2048 x 1536. Widefield 10x /20mm eyepieces. Binocular viewing head with 30° inclined eyepieces. 4x, 10x, 40xR & 100xR oil immersion objectives. 3 watt LED provides bright, white illumination for up to 50,000 hours

With Semi-plan Objectives:

#M10DB-S

With Plan Objectives:

#M10DB-P

Multi-Phase Objectives:

#M10DB-MP

IV Pump & Fluid Warmer



Mindray Benefusion VP3 Vet Infusion Pump:

- Can be used with generic IV sets
- Compact and light weight design for easy carrying and mounting
- Long runtime battery, for emergency and transfer
- Dynamic pressure display, prevents occlusion caused by small animal movement
- Sequential mode supports secondary infusion (piggyback)

#8610B-CTO-S01



Fluid Warmer with
LCD Display

#FW-300VET



Autoclaves & Ultrasonic Cleaners



EZ9 Automatic Autoclave 9"

[#EZ9](#)

EZ10 Automatic Autoclave 10"

[#EZ10](#)



EZ9Plus Automatic Autoclave 9"

[#EZ9PLUS](#)

EZ11Plus Automatic Autoclave 11"

[#EZ11PLUS](#)



Fully Automatic Vet Autoclave 10"

[#TVET 10E](#)

Fully Automatic Vet Autoclave 11"

[#TVET 11E](#)



Manual Autoclave 9"

[#2340M](#)

Manual Autoclave 10"

[#2540M](#)



M9 UltraClave Automatic Sterilizer

[#M9-043](#)

M11 UltraClave Automatic Sterilizer

[#M11-043](#)



UltraFast Automatic Sterilizer,
7.6" x 1.6"

[#M3-001](#)



Masteri 18L, Class B Steam Autoclave, 9.5" x 13.5"

[#AU2028V](#)

Masteri 23L, Class B Steam Autoclave, 9.5" x 18"

[#AU2033V](#)

Replacement filters

[#STE1028](#)



Ultrasonic Cleaner, 1 Gal

[#CSU1](#)

Ultrasonic Cleaner, 1 Gal
with Heater

[#CSU1H](#)

Ultrasonic Cleaner, 3 Gal

[#CSU3](#)

Ultrasonic Cleaner, 3 Gal
with Heater

[#CSU3H](#)



QuickClean:

Ultrasonic Cleaner 1.2 Gal

[#QC1-01](#)

Ultrasonic Cleaner 3.3 Gal

[#QC3-01](#)

Ultrasonic Cleaner 6.6 Gal

[#QC6-01](#)



Chamber Brite
Autoclave Cleaner
Box of 10 packets

[#CB0010](#)



Clean & Simple
Ultrasonic/Enzymatic
Tablets (box/64)

[#CS0064](#)

Stainless Steel, Cabinets and Caging



V-Top Surgery Table with Adjustable Hydraulic Column, 50"

[#100-4041-00](#)

V-Top Surgery Table with Adjustable Hydraulic Column, 60"

[#100-4071-00](#)

Flat-Top tables also available

IV Rod for Surgical Tables

[#200-3041-09](#)



Pedestal Exam Table

[#101-9640-00](#)



Mobile Electric Lift Table

[#103-0170-16](#)

With Midmark Platform Scale

[#110-3410-08](#)



Aluminum IV Pole,, 5 leg, 2 hook

[#52022](#)

Aluminum IV Pole, 5 leg, 4 hook

[#52023](#)



Mayo Instrument Table, Double Post

[#52027](#)



Faucet Mount Eye Wash Station

[#107-8492-00](#)



Medical Tray

[#52029](#)



Cage Bank 5 Unit Assembly Option 3


[#160-0105-20](#)



Hinge Mounted IV Rod



O₂ Therapy Door



Plastisol Coated Cage Floors

This is just a small selection of our products

Call for more options

Imaging Accessories



Imaging Chamber:

3 Mouse: [#078920341](#)

5 Mouse: [#078919969](#)



EquaFlow 3XL
Mouse/Rat:

[#EQF3001](#)



EquaFlow 5 Station
Mouse:

[#EQF5001](#)



EquaFlow 10 Station
Mouse:

[#EQF10001](#)



Optical Imaging
Chamber:

[#941440](#)



CX-IM Non-Rebreathing
System:

[#921439](#)

Rat Phantom:

Radiopaque (X-ray/CT) and Fluorescent
(Green, Red, and Cy5 Filter sets)

Radiopaque (X-ray/CT) and Luminescent
(Blue emission)

Luminescent (blue emission) and Fluorescent
(Green, Red, and Cy5 Filter Sets)



Euthanasia Supplies



Easy to use Preset Flow System (PFS) is designed to eliminate the variables the lead to non-compliance. Delivers accurately measured volumes for single housing cages, allowing any lab or facility to have consistent, reliable CO₂ delivery.

Preset Flow System for Mice and
Rats:

With Nut & Stem For Large Tanks:

[#934105](#)

With Yoke for Small Tanks:

[#934115](#)

For Wall Outlet:

[#934130](#)

CO₂ Flowmeter 0-12 LPM with DISS
Fittings:

[#FM-803C](#)

CO₂ Regulator for small tanks with 0-12
LPM flow meter:

[#M1-940-12FM](#)

CO₂ Regulator for large tanks with 0-12
LPM flow meter:

[#M1-320-12FM](#)



Universal Lid, Rat

8.5" x 9.75"

[#941445](#)



Universal Lid, Mouse

11.5" x 7.25"

[#941442](#)

Consumables



VaporGuard - 8 pk
Charcoal Filter

#931401



Amsorb Plus
Twin Pack

#AMAB3000



Somni W.A.G.
4 pack

#WG-15003



F/Air Canister
8 pack

#20522



WAG Canister
4 pack

#07-8909457



Charcoal
Prepack

#20510



VetLube
30mL Spray Bottle

#D14002



VaporGuard XL-R Refillable
Charcoal Filter (2 pack)

#931402

Bulk Carbon Container

#101500



Powered Autoclave Cleaner
Box of 10 packets

#CB0010



Redux Gel

#PL-200



Sting-Free Adhesive
Remover Spray and Wipes
for the fast, painless and
atraumatic removal of
adhesives.

50mL - **#TV-101**

200mL - **#TV-201**



StopLik wide 3" x 12" (12 / box)

#57605

StopLik narrow 1.62" x 12" (12 / box)

#57615

StopLik mixed box (6 wide, 6 narrow)

#57616



Clean & Simple - Ultrasonic/
Enzymatic Tablets (box/64)

#CS0064



Waste Gas Monitoring
Badge and Analysis Report

#931443



Sting-Free Barrier Spray
A durable, breathable spray
to protect the skin from
body, waste, enzymes and
adhesives

50mL - **#TV-104**

200mL - **#TV-204**



Monorapid Hand Sanitizer
1L with Pump

#AS-MAN-1000P

120mL Bottle

#AS-MAN-120



Smart-Y for Centesis
Package of 5

#AH7073



Disposable
Table Covers

#V501



6" x 30" 12 mil Latex
Sheet for rodent masks

#07-8917456

Examination gloves available in Latex, Vinyl or
Nitril, powder or powder free in sizes x-small to
x-large.

Bouffant caps 21" in a variety of colors and 24"
white

Call for Details



Reusable/Autoclavable Stapler Handle

#S1004

35 Wide Stapler Cartridge

#S1003-12





Single Pack

#SS101-1



SutureSeal™ is a long lasting, transparent bandage replacement made of a non-toxic biodegradable hydrogel. This is a specially formulated hydrogel designed for incision and wound protection in small and large animals. SutureSeal™ acts by forming a thin coating that protects the suture site or wound. It is Initially applied as a liquid and transforms into a flexible gel in about 90 seconds. SutureSeal™ forms a protective layer that prevents dirt and debris from entering and lasts for up to two weeks. SutureSeal™ is sterile, colorless, odorless and designed to enhance the healing process. The wound will heal and the polymer will dissolve naturally after two weeks.



AQUACEL® Hydrofiber® Wound Dressing, Sterile

AQUACEL® Hydrofiber® Dressing is a white, fibrous dressing derived from 100% pure sodium carboxymethylcellulose (NaCMC). Combines the look, feel and handling properties of gauze and leading alginates, but retains more fluid than alginates and almost three times more fluid than gauze when used under compression.

-Locks in wound exudate¹¹ and traps bacteria^{8-9, †}

-Micro-contours to the wound bed^{13, †}

-Responds to wound conditions by forming a cohesive gel

AQUACEL® EXTRA™ Hydrofiber® Wound Dressing

Same gelling benefits of the original AQUACEL® dressing, but is 9X stronger¹⁹ and 39% more absorbent¹⁹

AQUACEL® Foam Hydrofiber® Dressing, Sterile

A Soft absorbent foam pad, layered with an AQUACEL® interface that gels on contact, topped with a Waterproof/Bacteria Barrier that evaporates excess moisture.

AQUACEL® Surgical Hydrofiber® Cover Dressing

Combines Hydrofiber® Technology and hydrocolloid technologies as a single-product solution to the post-operative management of surgical incisions. Designed to enhance clinical outcomes and facilitate patient mobility.

KALTOSTAT® wound dressings, Sterile

A calcium-sodium alginate fiber dressing that can be used in moderately to heavily exuding wounds. Designed to absorb exudate and to help control minor bleeding.

Multiple sizes available. Most AQUACEL® dressings available with added silver (Ag) to provide rapid and sustained antimicrobial activity^{7, 9, 16, †}



7. Jones SA, Bowler PG, Walker M, Parsons D. Controlling wound bioburden with a novel silver-containing Hydrofiber® dressing. Wound Rep Reg. 2004;12:288-294.

8. Waring MJ, Parsons D. Physico-chemical characterisation of carboxymethylated spun cellulose fibres. Biomaterials. 2001;22:903-912.

9. Walker M, Hobot JA, Newman GR, Bowler PG. Scanning electron microscopic examination of bacterial immobilisation in a carboxymethylcellulose (Aquacel) and alginate dressings. Biomaterials. 2003;24(5):883-890.

13. Jones S, Bowler PG, Walker M. Antimicrobial activity of silver-containing dressings is influenced by dressing conformability with a wound surface. Wounds. 2005;17(9):263-270.

16. Bowler PG, Jones SA, Walker M, Parsons D. Microbicidal properties of a silver-containing Hydrofiber® dressing against a variety of burn wound pathogens. J Burn Care Rehabil. 2004;25:192-196.

19. Assessment of the Physical Properties of AQUACEL® EXTRA & AQUACEL Dressings. Scientific Background Report WHRI3461 TA214. 2011 Data on File, ConvaTec.

Accelerated Hydrogen Peroxide (AHP) is a synergistic blend of commonly used, safe ingredients that when combined with low levels of hydrogen peroxide produce exceptional potency as a germicide and performance as a cleaner. The benefits and efficacy of AHP have been validated by third party clinical studies conducted by scientific organizations and third party researchers that are recognized by government regulatory agencies in Canada and the U.S.



Prevail RTU(Ready to Use):

1 L

[#12101-001R](#)

4 L

[#12105-004R](#)



PREVail Wipes, 6"x7"

160 wipes

[#12221-6x7W](#)



Prevail Concentrate:

4 L

[#12305-004C](#)

20 L

[#12309-020C](#)

200 L

[#12311-200C](#)



Mixing Station:

1 Way - [A8031A-G](#)

2 Way - [A8032A-G](#)

3 Way - [A8033A-G](#)



Pump-Up Foamer:

1.5 Litre

[#S2100](#)

8 Litre

[#S2500](#)



Cordless Handheld
Electrostatic Sprayer

[#V2200](#)



Cordless Backpack
Electrostatic Sprayer

[#V3300](#)



4L hand dispensing pump

[#A8050](#)



Foaming Trigger Sprayer

[#A8019A](#)



PREempt HLD5 Instrument
High-Level Disinfectant

[#11505-004L](#)



PREempt CS20
Instrument Disinfectant

[#11405-004L](#)



No Drip Pump
For 200L Drum

[#A4600BKV](#)



Backpack Sprayer

[#S4100](#)

Battery Operated (Blue)

[#S4500](#)



Pure Oxygen Ultra
Shampoo Concentrate

[#A8334](#)



PREempt 1L
Sporocidal Liquid

[#11701-001L](#)



50 Test Strips

[#AHP1750](#)

Kennel Foamer
With 32oz Bottle

[#A8040B-32](#)

Cleaning Accessories



Long Handled Scrub Brush

[#F2760-XX](#)



Vikan Dustpan

[#F2820-XX](#)



20" Single Blade
Ultra Hygiene Squeegee

[#F2400-XX](#)



20" Foam Blade
Squeegee

[#F2430-XX](#)



1-3 Tool Wall Bracket

[#F2970-XX](#)

(4-6 Tool also available)



Ergonomic Scrub
Brush

[#F2740-XX](#)



Extra Large
Scrub Brush

[#F2730-XX](#)



Wall/Floor Scrub
Brush

[#F2870-XX](#)



51" Ultra Hygiene Handle

[#F2640-XX](#)

(59" also available)



Detail Brush

[#F2810-XX](#)



6" Mini Handle

[#F2610-XX](#)



Narrow Utility Brush

[#F2780-XX](#)



High-Low Brush

[#F2860-XX](#)

Your Local Benson Rep.

BC & Southern Alberta:

Donna Locke, BSc., RVT, RLAT

Clinical Product Specialist

1-800-563-3859 ext 272

dlocke@bensonmedical.ca

Northern Alberta, Saskatchewan & Manitoba:

Jodi Walchuk, RVT

Clinical Product Specialist

1-800-563-3859 ext 273

jwalchuk@bensonmedical.ca

Ontario:

Krista Holdsworth, AHT, RLAT

Clinical Product Specialist

1-800-563-3859 ext 275

kholdsworth@bensonmedical.ca

Ottawa & Western Quebec:

Mark Landry, AHT

Clinical Product Specialist

1-800-563-3859 ext 281

mlandry@bensonmedical.ca

Atlantic Provinces, Eastern Quebec:

Shauna Cormier, RVT

Clinical Product Specialist

1-800-563-3859 ext 282

scormier@bensonmedical.ca

ANAESTHESIA VAPORIZER SERVICE

Manufacturer's comprehensive overhaul and calibration includes:



- * Complete disassembly and cleaning
- * Internal and external inspection
- * Replacement of internal wicks and seals
- * Calibration and adjustment of Thermo-Control
- * Leak test and performance verification
- * Rental vaporizers provided upon request

On site vaporizer output verification is **NOT** a substitute for service & calibration

Vaporizers ARE NOT on site serviceable



TEC III



TEC IV



TEC V



OHIO

We also do Annual Preventative Anaesthetic Machine Maintenance

Veterinary service companies are highly unregulated, it is important to be aware that Benson Medical has been servicing anaesthetic equipment since 1974 and maintains the following important business and service practices;

- * We maintain a Health Canada establishment license.
- * We are authorized Canadian distributor for your most prominent anaesthesia equipment.
- * We maintain liability insurance.
- * We offer a comprehensive service report.



benson
medical industries inc.

www.bensonmedical.ca
1.800.563.3859